


UHPLC Column for Hydrophilic Compound Separation

COSMOSIL 1.8PBr

COSMOSIL PBr features a pentabromobenzyl functional group, which is suited to the retention of hydrophilic compounds in reversed-phase mode. This product is an ultra-high-performance liquid chromatography (UHPLC) with a particle size of 1.8 μm , enabling rapid separation of hydrophilic compounds under familiar reversed-phase analytical conditions.


POINT

- 01 Retain hydrophilic compounds with reversed-phase chromatography**
- 02 Significantly improved separation compared to ODS columns**
- 03 High-speed separation with UHPLC (1.8 μm particles)**



Stationary Phase Structure

Ultra-high performance liquid chromatography (UHPLC) is a technology for speeding up analysis by using columns packed with smaller silica gel particles than standard HPLC and flowing mobile phase at a high rate. Analyses by UHPLC are increasing in number every year, and it is sure to continue to be an important technology in years to come.




Number of Publications Using UHPLC (2005-2024)
Number of publications with 'UPLC' or 'UHPLC' in the title, abstract, or keywords (searched via <https://www.sciencedirect.com/>)

Product description

Rapid separation of hydrophilic compounds using UHPLC

Rapid separation of purine metabolites


Using COSMOSIL PBr, the analysis of 9 purine nucleoside metabolites, including uric acid and allantoin, required approximately 40 minutes. Changing to COSMOSIL 1.8PBr, packed with smaller particles, allowed the analysis to be completed in approximately 10 minutes with equivalent separation.




Product description (cont.)

Hydrophilic compound retention by interaction with bromine


COSMOSIL PBr columns can be used in the same reversed-phase mode as ODS columns; however, due to the high polarizability of bromine atoms, it exhibits unique retention and separation behavior. This is attributed to dispersion interactions and the recently recognized halogen bonding, which differ from conventional ODS columns. These interactions occur between highly polarizable atoms and π -electrons, as well as oxygen or nitrogen atoms possessing lone electron pairs. As a result, PBr columns can retain hydrophilic compounds that are difficult to retain on ODS columns, and enable the analysis of difficult compounds.



Imidazole dipeptide analysis



Nicotinamide metabolite analysis




◆ Differences in separation characteristics due to stereoselectivity compared to ODS columns

With COSMOSIL PBr columns, selectivity derives not only from π -interactions of aromatic rings and dispersion interactions with bromine, but also structural configuration of the analytes. These effects can lead to improved separation on PBr columns, even when retention on ODS columns is sufficient.

* Taniguchi A, et al. Chromatography. 2025;46:55.


● Malodorous aldehyde derivative analysis

In the analysis of odorous aldehyde derivatives regulated under Japan's Offensive Odor Control Law, separation of certain isomers can be challenging. When using conventional reversed-phase columns, complex analytical conditions, such as a three-component mobile-phase gradient, have been used. In contrast, analysis can be performed using 1.8PBr under simpler isocratic conditions using a two-component mobile phase.



◆ Differences in separation characteristics compared to PFP columns

PBr columns and pentafluorophenyl (PFP) columns both contain halogenated benzene. However, because PFP relies on dipole–dipole interactions, while PBr retains compounds primarily through dispersion interactions, their separation characteristics differ significantly. Both PFP and PBr columns showed improved separation compared to ODS columns, but the elution order was different. PFP exhibited strong retention of the fluorinated groups, $-\text{CF}_3$ and $-\text{SF}_5$, while PBr exhibited strong retention of the heavy halogen groups, $-\text{Br}$ and $-\text{I}$, through dispersion interactions.




◆ Material properties


Packing Material	1.8PBr		
Silica Gel	Fully Porous Spherical High-Purity Silica	Surface Area (m^2/g)	340
Particle Size (μm)	1.8	Bonded Phase	Pentabromobenzyl group
Pore Size (nm)	12	Bonding Type	Monomeric

Applications


● Purity of synthetic peptides




● Separation of isomers of highly hydrophilic compounds



● Separation of adenosine phosphates



● Separation of sugars



Related resources

For additional product information, please visit our website. You can access it using the QR codes below.

COSMOSIL PBr Series



Information on other UHPLC series




COSMOSIL Application
(Search for analysis examples)



Note on connector type:

Our UHPLC columns use the same connectors as Waters UPLC (UHPLC) columns. This is different from our conventional COSMOSIL columns, which use the conventional Waters HPLC-compatible connectors. Attempting to connect an unsuitable fitting may result in it becoming stuck in the column.



All product names, trademarks, and registered trademarks are the property of their respective owners. Use of these names does not imply any affiliation or endorsement.

Ordering information

COSMOSIL 1.8PBr Packed Column

Column Size I.D. x Length (mm)	Product Number	Column Size I.D. x Length (mm)	Product Number	Column Size I.D. x Length (mm)	Product Number
2.1 x 30	22739-41	2.1 x 100	22742-81	3.0 x 50	22744-61
2.1 x 50	22740-01	2.1 x 150	22743-71	3.0 x 100	22745-51
2.1 x 75	22741-91			3.0 x 150	22746-41

COSMOSIL is a registered trademark of Nacalai Tesque, Inc.

For research use only, not intended for diagnostic or drug use.