



COSMOSIL

Fatty Acid Analysis by HPLC

Technical Note

Current interest in fatty acids and trans-fatty acids in food materials has been increasing. Fatty acid is an organic acid found in fats, oils, and lipids. It can be classified according to chain length, short chain (2-4 carbons), middle chain (5-10 carbons), and long chain (more than 11 carbons). Although gas chromatography is the predominant technique used for fatty acid analysis, high-performance liquid chromatography (HPLC) plays an important role in applications such as the geometrical isomer separation. By using both HPLC and GC, a better fatty acid profile can be obtained.



		High-Performance Liquid Chromatography (HPLC)		Gas Chromatography (GC)
		COSMOSIL C ₁₈	COSMOSIL Cholester	
Free fatty acids	Short fatty acids	Good	(HILIC)	Poor
	Middle and long fatty acids	Excellent		Poor
Fatty acid ester	Short fatty acids	Fair		Poor
	Middle and long fatty acids	Good		Excellent
Position or geometrical isomers of unsaturated fatty acids		Fair	Good	Excellent
High sensitive analysis		Fair [labeling] Excellent		Excellent
Preparative separation		Excellent		Poor
Lipids composed of complex fatty acid mixture		Using both HPLC and GC to achieve better fatty acid profiling (*1)		

*1: Wakako Tsuzuki and Kaori Ushida, "Preparative Separation of *cis*- and *trans*- Isomers of Unsaturated Fatty Acid Methyl Esters Contained in Edible Oils by Reversed-Phase High-Performance Liquid Chromatography" *Lipids* (2009) 44: 373-379

1. Unsaturated Fatty Acid Analyses by HPLC

(1) Geometrical Isomer Separations (*cis/trans*)

Geometrical isomers (*cis/trans*) are difficult to separate by a C₁₈ column due to their similar hydrophobicity. COSMOSIL Cholester column can achieve better separation due to higher molecular-shape selectivity.

• 18 :1 *cis*-9 / *trans*-9 Fatty Acid

COSMOSIL Application Data

Column: 4.6mm I.D.-150mm
 Mobile phase: 0.05%TFA-90%Methanol
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD

Sample: 1; Oleic Acid (3.0 μg)
 2; Elaidic Acid (3.0 μg)
 3; Stearic Acid (3.0 μg)

5C₁₈-MS-II

Cholester



NACALAI TESQUE, INC
AP-1038

• 18 :1 *cis*-11 / *trans*-11 Fatty Acid

COSMOSIL Application Data

Column: 5C₁₈-MS-II
 Column size: 4.6mm I.D.-150mm
 Mobile phase: 0.05%TFA-90%Methanol
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD

Sample: 1; *cis*-Vaccenic Acid (3.0 μg)
 2; *trans*-Vaccenic Acid (3.0 μg)



NACALAI TESQUE, INC

AP-1039

(2) Positional Isomers


COSMOSIL Cholester column with high molecular-shape selectivity offers improved separation for positional isomers. Labeling can increase detection sensitivity.

• Labeled 18 : 3 (n-3) / 18 : 3 (n-6) Fatty Acid Ester

COSMOSIL Application Data

Column: 5C₁₈-MS-II
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/ H₂O = 90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254nm

Sample: 1; Linolenic Acid
p-Bromophenacyl Ester
 2; γ-Linolenic Acid
p-Bromophenacyl Ester




NACALAI TESQUE, INC

AP-0209

COSMOSIL Application Data

Column: Cholester
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Methanol/ H₂O = 90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254nm

Sample: 1; Linolenic Acid *p*-Bromophenacyl Ester
 2; γ-Linolenic Acid *p*-Bromophenacyl Ester



NACALAI TESQUE, INC

AP-0211

2. Free Fatty Acid Analyses by HPLC

Methyl esterification of free fatty acids is required for GC analyses to obtain better peak shapes and lowered boiling points. Derivatization is not required for HPLC analyses.



(1) Short-chain Free Fatty Acids

Short-chain free fatty acids have less retention on C₁₈ columns. Therefore COSMOSIL C₁₈-PAQ column (100% water compatible) and COSMOSIL HILIC column are better for short-chain free fatty acids analyses.

• Propanoic Acid (3:0)

COSMOSIL Chromatogram Index

Sample: Propionic Acid
 CAS No.: [79-09-4]
 Molecular formula: C₃H₆O₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0μl
 Retention time: 14.85min
 Capacity factor: 4.24





NACALAI TESQUE, INC

• Acrylic Acid (3:1)

COSMOSIL Chromatogram Index

Sample: Acrylic Acid
 CAS No.: [79-10-7]
 Molecular formula: C₃H₄O₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0μl
 Retention time: 15.05min
 Capacity factor: 4.28




NACALAI TESQUE, INC

• Butyric Acid (4:0)

COSMOSIL Chromatogram Index


Sample: *n*-Butyric Acid
 CAS No.: [107-92-6]
 Molecular formula: C₄H₈O₂
 Column: 5C18-MS-II
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/20mmol/l Phosphate Buffer (pH2.5) =40/60
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4x10⁻⁵ RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 3.79min
 Capacity factor: 1.24



• Crotonic Acid (4:1)

COSMOSIL Chromatogram Index

Sample: Crotonic Acid
 CAS No.: [3724-65-0]
 Molecular formula: C₄H₆O₂
 Column: 5C18-MS-II
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/20mmol/l Phosphate Buffer (pH2.5) =30/70
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4x10⁻⁵ RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 4.21min
 Capacity factor: 1.46




(2) Middle or Long-chain Free Fatty Acids

COSMOSIL C₁₈-MS-II column was used for non-derivatized middle or long-chain free fatty acid analyses.

• Capronic Acid (6:0)

COSMOSIL Chromatogram Index


Sample: *n*-Capronic Acid
 CAS No.: [142-62-1]
 Molecular formula: C₆H₁₂O₂
 Column: 5C18-MS-II
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/20mmol/l Phosphate Buffer (pH2.5) =60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4x10⁻⁵ RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 4.63min
 Capacity factor: 1.72



• Pelargonic Acid (9:0)

COSMOSIL Chromatogram Index


Sample: Pelargonic Acid
 CAS No.: [112-05-0]
 Molecular formula: C₉H₁₈O₂
 Column: 5C18-MS-II
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/20mmol/l Phosphate Buffer (pH2.5) =80/20
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4x10⁻⁵ RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 4.29min
 Capacity factor: 1.56



• Lauric Acid (12:0)

COSMOSIL Chromatogram Index


Sample: Lauric Acid
 CAS No.: [143-07-7]
 Molecular formula: C₁₂H₂₄O₂
 Column: 5C18-MS-II
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/20mmol/l Phosphate Buffer (pH2.5) =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4x10⁻⁵ RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 4.18min
 Capacity factor: 1.51








• Pentadecanoic Acid (15:0)

COSMOSIL Chromatogram Index

Sample: *n*-Pentadecanoic Acid
 CAS No.: [1002-84-2]
 Molecular formula: C₁₅H₃₀O₂
 Column: 5C18-MS-II
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/20mmol/l Phosphate Buffer (pH2.5) =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4x10⁻⁵ RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 2.5µl
 Retention time: 7.97min
 Capacity factor: 3.79



Specifications

Packing Material	C ₁₈ -MS-II	C ₁₈ -PAQ	Cholester	πNAP	PYE	HILIC
Bonded Phase Structure						
Bonded Phase	Octadecyl Group	Octadecyl Group	Cholesteryl Group	Naphtylethyl Group	Pyrenylethyl Group	Triazole
Main Interaction	Hydrophobic Interaction	Hydrophobic Interaction	Hydrophobic Interaction Molecular Shape Selectivity	Hydrophobic Interaction π-π Interaction	Hydrophobic Interaction π-π Interaction Stereoselectivity Charge-transfer Interaction	Hydrophilic Interaction
Features	<ul style="list-style-type: none"> Multi-purpose C₁₈ column for separation of the widest range of compounds 	<ul style="list-style-type: none"> Compatible with 100% water based mobile phase. Suitable for hydrophilic compounds. 	<ul style="list-style-type: none"> Specialty for structural isomers Usable under the same condition as C₁₈. 	<ul style="list-style-type: none"> Stronger π-π interaction than phenyl column 	<ul style="list-style-type: none"> The most powerful π-π interaction 	<ul style="list-style-type: none"> Suitable for highly-polar compounds Ion-pair reagent is not required
Product Code Column Size: 4.6 x 150	38019-81	02486-71	05976-61	08085-41	37837-91	07056-51

COSMOSIL Application

COSMOSIL Application has more than 7,000 applications using COSMOSIL columns. Setting optimal HPLC experimental parameters is the one of the most important processes that requires experience and time. COSMOSIL Application provides you with sample analysis conditions with widely used ODS columns and other specialty columns.

Visit our web site at <http://www.nacalai.co.jp/cosmosil/data/csmosrctop.cfm?lc=E> or type "COSMOSIL Application".

Only show new applications added on Dec. 20th 2011

Category:

Column name:

Sample Name: begins with

CAS number: (ex:498-02-2)

Particle Size:

Result/Page:

Applications are search by

1. Sample Category
2. Sample Name
3. CAS No.,
4. Column Name
5. Particle Size

Search Result

COSMOSIL Application


Search condition [Column:πNAP]

[TOP]

Results 24 (1-20) [Next]

Data No.	Data Name	Sample	Particle Size (μm)	Column	CAS No.
AP-1206	Dichlorophenol	2,3-Dichlorophenol	5	πNAP	576-24-9
		2,4-Dichlorophenol			120-83-2
		2,5-Dichlorophenol			583-78-8
		2,6-Dichlorophenol			87-65-0

COSMOSIL Application



For research use only, not intended for diagnostic or drug use.

nacalai tesque
The quality for certainty.

NACALAI TESQUE, INC.
Nijo Karasuma, Nakagyo-ku, Kyoto 604-0855 JAPAN
TEL : +81-(0)75-251-1730
FAX : +81-(0)75-251-1763
Website : www.nacalai.com
E-mail : info.intl@nacalai.com