


Data courtesy of Prof. Ryuji Kato, Ph D.,  
Laboratory of Cardiovascular Pharmacotherapy and Toxicology,  
Osaka University of Pharmaceutical Sciences

COSMOSIL πNAP uses a naphthalene-based bonded phase, and exhibits different retention behavior from C<sub>18</sub>. In this note, πNAP is used to analyze cisplatin (CDDP) and its mono-chloro, none-chloro (hydration complex), and OH-dimer forms in water.

Standard C<sub>18</sub> columns have difficulty separating compounds of similar hydrophobicity. Performance may be improved by using Cholester, a column with improved molecular shape recognition.

### (1) Experiment Summary

Cisplatin (CDDP) in water converts to its mono-chloro, hydration complex, and OH-dimer forms, and it is thought that they exist in an equilibrium state (fig. 1). Each of these forms was separated and quantified.



**Fig. 1. Equilibrium reactions of cisplatin (CDDP) in water**  
In water, CDDP loses  $\text{Cl}^-$ , which is replaced with water. Depending on the pH, it may be replaced by a hydroxide ion. The mono-chloro form may lose another  $\text{Cl}^-$ , and, in the same way, the water molecules (one or both) may be replaced with a hydroxide ion. These forms are in equilibrium.

### (2) Quantification of Cisplatin (CDDP) and its Mono-Chloro, None-Chloro, and OH-Dimer Forms

- Comparing a C<sub>18</sub> column, COSMOSIL Cholester, and πNAP under identical conditions





When measuring the OH-dimer nitrate form, the dissociated  $\text{NO}_3^-$  peak elutes before the OH-dimer peak.

C<sub>18</sub>-MS-II could not separate the two peaks, and Cholester achieved partial separation, but could not separate to baseline. πNAP, however, was able to separate the compounds completely.

#### COSMOSIL Application Data

Column: COSMOSIL \*\*  
Column size: 4.6 mm I.D. × 250 mm  
Mobile phase: 100 mmol/l  $\text{NaClO}_4$ /Acetonitrile/100 mmol/l  $\text{HClO}_4$   
= 290/10/3  
Flow rate: 1.0 ml/min  
Temperature: 37 °C  
Detection: UV225 nm


Sample: Cisplatin OH-dimer nitrate (2 mmol/l)  
1; Nitrate  
2; Cisplatin OH-dimer  
Injection Vol. 20 μl



Data courtesy of Ryuji Kato, Ph D., Laboratory of Cardiovascular Pharmacotherapy and Toxicology,  
Osaka University of Pharmaceutical Sciences

## ■ Analysis Using COSMOSIL πNAP

Below is a typical chromatogram; the peak separation was excellent. The calibration curve for cisplatin (CDDP) was linear in the range of 0.01 to 4 mM, and the CV values for within-day and between-day variation were within 5%.



### (3)Observations

By using COSMOSIL πNAP, it is possible to separate and quantify cisplatin (CDDP) and its mono-chloro, none-chloro, and OH-dimer forms using HPLC.

### (4)References

Kato R. et al. A Novel Analytical Method of Cisplatin Using the HPLC with a Naphthylethyl Group Bonded with Silica Gel (πNAP) Column. *Biol Pharm Bull.* **40**(3), 290-296(2017).

For research use only, not intended for diagnostic or drug use.