

Fine-powdered Agarose

For Agarose Electrophoresis







One of the advantages of the fine-powdered agarose is the increased solubility. Not only easy to handle and weigh, it saves your precious time. The most important feature is the high resolution of the gel it produces.

Features

High soluble: smaller average particle size for easy dissolution
Simple: easy weighing operation
Clear: sharp and clean electrophoresis result


Solubility comparison - Particle size and solubility

Solubility tests have been done using microwave oven (400 W).
Photo image: microscope (x25)



Agarose:	Conventional	Fine-powdered	A Company	B Company
Solubility (speed):	Fast	Very fast	Very fast	Normal
Time:	2 min 50 sec	1 min 53 sec	1 min 55 sec	3 min 10 sec

Image comparison



Agarose:	Conventional	Fine-powdered	A Company	B Company
Electrophoresis conditions :				
Gel	: 1.5% Agarose	Sample	① Perfect DNA Markers, 0.5-12kbp	0.1 μg/μl
Buffer	: 1xTAE Buffer		② λ/HindIII digest φ x174/HincII digest	0.1 μg/μl
Time	: 75 min		③ 100bp DNA Ladder	0.12 μg/μl
Dyeing method	: EtBr			

Specification

Type :	≥ 1kbp
Sulfate (%) :	≤ 0.2
Gel Strength :	≥ 2,500 g/cm ² (at 1.5%)
Gel Point (°C) :	36 ± 1.5
Electroendosmosis (-mr) :	0.09-0.13

Ordering information

Product name	Grade	Storage	Code No.	PKG Size
Fine-powdered Agarose	SP	RT	02468-66	100 g
			02468-95	500 g

Ordering Information

Product name	Grade	Storage	Code No.	PKG Size
Agarose for ≥ 1 kbp fragment	SP	RT	01145-32	25 g
Type : ≥ 1 kbp			01145-74	100 g
Sulfate (%) : ≤ 0.5				
Gel Strength : $\geq 2,000$ g/cm ² (at 1.5%)				
Gel Point (°C) : 35 - 37				
Electroendosmosis (-mr) : ≤ 0.12				
Agarose for ≥ 1 kbp fragment	SP	RT	01163-92	25 g
Type : ≥ 1 kbp			01163-76	100 g
Sulfate (%) : ≤ 0.2			01163-05	500 g
Gel Strength : $\geq 2,500$ g/cm ² (at 1.5%)				
Gel Point (°C) : 36 ± 1.5				
Electroendosmosis (-mr) : 0.09-0.13				
Agarose-RE for ≥ 1 kbp fragment, for Restriction and Ligation	SP	RT	01149-92	25 g
Type : ≥ 1 kbp			01149-76	100 g
Sulfate (%) : ≤ 0.2			01149-05	500 g
Gel Strength : $\geq 2,500$ g/cm ² (at 1.5%)				
Gel Point (°C) : 36 ± 1.5				
Electroendosmosis (-mr) : 0.09 - 0.13				
DNase, RNase tested				
Agarose for 150-1,500bp fragment	SP	RT	01153-22	25 g
Type : 150 - 1,500bp			01153-64	100 g
Sulfate (%) : ≤ 0.1				
Gel Strength : $\geq 2,000$ g/cm ² (at 1.5%)				
Gel Point (°C) : ≤ 36.5				
Electroendosmosis (-mr) : ≤ 0.12				
DNase, RNase tested				
Agarose for 50-800bp fragment	SP	RT	01147-12	25 g
Type : 50 - 800bp			01147-96	100 g
Sulfate (%) : ≤ 0.1				
Gel Strength : ≥ 750 g/cm ² (at 1.5%)				
Gel Point (°C) : 30				
Electroendosmosis (-mr) : ≤ 0.12				
DNase, RNase tested				
Agarose-LE, Classic type	SP	RT	01157-82	25 g
Type : LE, Classic			01157-66	100 g
Sulfate (%) : ≤ 0.2			01157-95	500 g
Gel Strength : $\geq 2,500$ g/cm ² (at 1.5%)				
Gel Point (°C) : 36 ± 1.5				
Electroendosmosis (-mr) : 0.09 - 0.13				
Agarose-ME, Classic type	SP	RT	01158-72	25 g
Type : ME, Classic			01158-56	100 g
Sulfate (%) : ≤ 0.25			01158-85	500 g
Gel Strength : $\geq 2,000$ g/cm ² (at 1.5%)				
Gel Point (°C) : 36 ± 1.5				
Electroendosmosis (-mr) : 0.16 - 0.19				

For research use only, not intended for diagnostic or drug use.

NACALAI TESQUE, INC.

Nijo Karasuma, Nakagyo-ku, Kyoto 604-0855, JAPAN

TEL: +81-75-251-1730

FAX: +81-75-251-1763

Web site: <http://www.nacalai.com>

E-mail: info.intl@nacalai.co.jp