

Protocol 101

Immobilization of Proteins on COOH Beads (EDC method)

Materials

1) Beads and Proteins

COOH beads (Product No.: TAS8848N1140): 1 mg

(Functional groups: 200–300 nmol/mg)

Proteins: 50 µg (= 1 nmol / 50 kDa of protein)

2) Reagents

2-Morpholinoethanesulfonic acid (MES)

2-[4-(2-Hydroxyethyl)-1-piperazinyl]ethanesulfonic acid (HEPES)

Sodium hydroxide (NaOH)

Potassium chloride (KCl)

Ethylenediaminetetraacetic acid (EDTA)

Tween 20

1-Ethyl-3-(3-dimethylaminopropyl)carbodiimide hydrochloride (EDC·HCl) (M.W. 191.70)

2-Aminoethanol (Ethanolamine) (M.W. 61.08)

3) Apparatus

Micro centrifuge (HITACHI CF15RX2)

Microtube Mixer (TOMY MT-360)

Method

1) Preparation of buffers

- **Binding buffer**

25 mM MES-NaOH (pH 6.0)

- **Washing buffer**

10 mM HEPES-NaOH (pH 7.9)

50 mM KCl

1 mM EDTA

0.1 % Tween 20

- **Masking solution**

1 M Ethanolamine

10 mM HEPES-NaOH (pH 7.9)

50 mM KCl

1 mM EDTA

2) Immobilization of Proteins on COOH Beads

Start with 1 mg of COOH beads

↓

Wash beads with H₂O, 3 times

(Centrifuge at 15,000 rpm for 5 min at 4 °C and remove the supernatant)

↓

Add 200 µl of 20 mg/ml EDC·HCl solution in H₂O

↓

Mix for 30 min at 4 °C by using Microtube Mixer

↓

Centrifuge at 15,000 rpm for 5 min at 4 °C and remove the supernatant

↓

Add 500 µl of cold H₂O and mix rapidly

↓

Centrifuge at 15,000 rpm for 5 min at 4 °C and remove the supernatant

↓

Add 500 µl of cold Binding buffer and Mix rapidly

↓

Centrifuge at 15,000 rpm for 5 min at 4 °C and remove the supernatant

↓

Add 50 µg of proteins in 50 µl Binding buffer

↓

Mix for 2 h at 4 °C by using Microtube Mixer

↓

Add 5 µl of Masking solution

↓

Mix for 16-20 h (over night) at 4 °C by using Microtube Mixer

↓

Wash with 1 ml of Washing buffer, 5 times

(Centrifuge at 15,000 rpm for 5 min at 4 °C and remove the supernatant)

↓

Resuspend in Washing buffer, and store at 4 °C