


## Protocol for detecting target protein

---

### ELISA by immuno-plate

1. Add suitable antigen solution (i.e. SARS-CoV-2 Spike protein, conc. 0.1  $\mu\text{g/ml}$ ~5.0  $\mu\text{g/ml}$ ) into immuno-plate such as MaxiSorp™ and incubate for O/N at 4°C
2. Wash each well with maximum capacity of PBS-T for five times
3. Add blocking solution such as Blocking One (Nacalai, #03953-95) and incubate for 2 hrs at RT
4. Remove blocking solution and apply this product (1:400-1:4,000) or your antibody samples
5. Incubate sample for 1hr at RT
6. Repeat step 2
7. Add enzyme-conjugated secondary anti-human antibody solution and incubate for 1hr at RT
8. Repeat step 2
9. Add substrate solution such as TMB, and incubate for 5 to 30 min at RT
10. Stop reaction by acid solution such as phosphoric acid or sulfuric acid
11. Read suitable absorbance, and analyze data




#### ELISA

Standard curve of this product against SARS-CoV-2 Spike protein

### Western blot

1. Prepare ~100 ng of protein samples for PAGE
2. Perform PAGE including protein standard such as Protein ladder One+ triple color (Nacalai, #19593-25)
3. Transfer segregated proteins into PVDF membrane
4. Soak PVDF into blocking solution such as Bullet Blocking One for Western Blotting (Nacalai, #13779-56) for 5 min~ at RT
5. Apply this product (1:500) or your antibody samples on PVDF for 1 hr
6. Wash with TBS-T for three times
7. Apply enzyme-conjugated secondary anti-human antibody solution and incubate for 1 hr at RT
8. Repeat step 6
9. Visualize target proteins using Chemi-Luminescent substrate such as Chemi-Lumi One Super (Nacalai, #02230-14) or Chemi-Lumi one Ultra (Nacalai, #11644-24)
10. Detect signals and analyze data



#### Western blot

Western blot analyses of this product against SARS-CoV-2 S-protein (S) or N-protein (N)