JANUARY 2019

SUMMARY :

REVIEW

Immune Checkpoint Blockade: InestimAble advances

PRODUCTS

Immune Checkpoint Anti-human Antibodies

- Anti-hCTLA4 Isotype Family
- Anti-hPD1 Isotype Family
- Anti-hPD-L1 Isotype Family

ADCC Reporter Cell-Based Assay

- Jurkat-Lucia™ NFAT-CD16 Cells
- Raji-hCTLA4 Cells
- Raji-hPD-1 Cells
- Raji-hPD-L1 Cells

pTRIOZ Plasmids

Zeocin™

InvivoGen

ASIA

Unit 106, 1F,
15W Phase 3, Hong Kong Science Park,
Pak Shek Kok,
Hong Kong

T: +852 3622 3480
F: +852 3622 3483

info.hk@invivogen.com
www.invivogen.com

Immune Checkpoint Blockade: InestimAble Advances

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
1. Kavelansky J and Pavlick A.C. 2017. Beyond checkpoint inhibitors: the next generation of immunotherapy in oncology. *AJHO*. 13(2):9.
2. Ribas A. and Wolchok J.D. 2018. Cancer immunotherapy using checkpoint blockade. *Science*. 359:1350.
3. Wei, S.C. et al. 2018. Fundamental mechanisms of immune checkpoint blockade therapy. *Cancer Discov*. 8(9):1069.
4. Quast I. et al. 2017. Regulation of antibody effector functions through IgG Fc N-glycosylation. *Cell. Mol. Life. Sci.* 74(5):837.
5. Almagro, J.C. et al. 2018. Progress and challenges in the design and clinical development of antibodies for cancer therapy. *Front. Immunol.* 8:1751.
6. Whang, X. et al. 2018. IgG Fc engineering to modulate antibody effector functions. *Protein Cell*. 9:63.
7. Donini, C. et al. 2018. Next generation immune-checkpoints for cancer therapy. *J Thorac Dis.* 10 (suppl 13): S1581.
8. André P. et al. 2018. Anti-NKG2A mAb is a checkpoint inhibitor that promotes anti-tumor immunity by unleashing both T and NK cells. *Cell*. 175:1731.
9. Marshall H.T. and Djamgoz B.A. 2018. Immuno-oncology: emerging targets and combination therapies. *Front. Oncol.* 8:315.

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
PRODUCT

CAT. CODE

Anti-hCTLA4 antibodies, variable regions of ipilimumab

Anti-hCTLA4-hIgG1	hctla4-mab1
Anti-hCTLA4-hIgG1Nq	hctla4-mab12
Anti-hCTLA4-hIgG1fut	hctla4-mab13
Anti-hCTLA4-hIgG2	hctla4-mab2
Anti-hCTLA4-hIgG4 (S228P)	hctla4-mab14
Anti-hCTLA4-hIgA2	hctla4-mab7

Anti-hPD1 antibodies, variable regions of nivolumab

Anti-hPD1-Ni-hIgG1	hpd1ni-mab1
Anti-hPD1-Ni-hIgG1Nq	hpd1ni-mab12
Anti-hPD1-Ni-hIgG1fut	hpd1ni-mab13
Anti-hPD1-Ni-hIgG2	hpd1ni-mab2
Anti-hPD1-Ni-hIgG4 (S228P)	hpd1ni-mab114
Anti-hPD1-Ni-hIgA2	hpd1ni-mab7

Anti-hPD1 antibodies, variable regions of pembrolizumab

Anti-hPD1-Pem-hIgG1	hpd1pe-mab1
Anti-hPD1-Pem-hIgG1Nq	hpd1pe-mab12
Anti-hPD1-Pem-hIgG2	hpd1pe-mab2
Anti-hPD1-Pem-hIgG4 (S228P)	hpd1pe-mab14
Anti-hPD1-Pem-hIgA2	hpd1pe-mab7

Anti-hPD-L1 antibodies, variable regions of atezolizumab

Anti-hPD-L1-hIgG1	hpd1l-mab3
Anti-hPD-L1-hIgG1 (N298A)	hpd1l-mab12
Anti-hPD-L1-hIgG1fut	hpd1l-mab13
Anti-hPD-L1-hIgG2	hpd1l-mab2

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Immune Checkpoint Inhibitor Target Cells

- Raji-hCTLA4 Cells
- Raji-hPD-1 Cells
- Raji-hPD-L1 Cells

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

ADCC Reporter Cell-Based Assay


ADCC reporter cell-based assay procedure. 1- **ADCC Induction**: the target cells expressing surface antigens are incubated with specific mAbs for 1 hour before the addition of Jurkat-Lucia™ NFAT-CD16 effector cells at a determined target:effector ratio for 6 hours. 2- **Addition of Lucia luciferase detection reagent QUANTI-Luc™**: detection reagent is added to supernatant samples. 3- **ADCC measure**: bioluminescent signal is detected using a luminometer.

Example of a comparison of ADCC potency for native and engineered isotypes of an anti-human CTLA4 using Raji-hCTLA4 target cells. Raji-hCTLA4 cells were incubated with gradient concentrations of Anti-hCTLA4 or Anti-β-galactosidase mAbs for 1 hour. Jurkat-Lucia™ NFAT-CD16 effector cells were then co-incubated with targets cells for 6 hours. NFAT activation, reflecting the induced ADCC response, was assessed by determining Lucia luciferase activity in the supernatant using QUANTI-Luc™. Percentages of the maximal response normalized to the IgG1 isotype are shown. See left page for each isotype description.

● IgG1 ● IgG2
● IgG1fut ● IgG4(S228P)
● IgG1NQ ● IgA2 ● β-gal

Increased ADCC activity mediated by IgG1 compared to IgG1fut (non-fucosylated): Raji-hCTLA4, -hPD-1, and -hPD-L1 cells were incubated with Jurkat-Lucia™ NFAT-CD16 effector cells and corresponding IgG1 or IgG1fut specific mAbs. The data represent the EC₅₀ for each antibody.

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
- **Light chain cassette**: FerL promoter + MCS + constant (C) IgG L chain
- **Heavy chain cassette**: FerH promoter + MCS + constant (C) IgG H chain
- **Zeocin™ selection cassette**: enhancer + composite promoter + *Sh Ble* gene

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
pTRIOZ general plasmid map. pTRIOZ contains three distinct cassettes with composite promoters for the optimized expression of the H and L chains, as well as the *Sh ble* gene conferring resistance to Zeocin™. It contains two unique multiple cloning sites (MCS) for the insertion of the variable H and L chains.

Comparison between pTRIOZ and pFUSE mAb yield. For the production of anti-PDL-1, an optimized 1:1 ratio of the pFUSE heavy chain and pFUSE light chain plasmids was tested alongside pTRIOZ. The maximum yield is normalized to the pFUSE ratio.

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PRODUCT	QUANTITY	CAT. CODE
Zeocin™ (solution)	1g (10 x 1ml)	ant-zn-1
Zeocin™ (solution)	5g (50 x 1ml)	ant-zn-5
Zeocin™	Bulk	Available upon enquiry

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