MARCH 2018

SUMMARY :

REVIEW

Inflammasomes: connecting innate and adaptive immunity

PRODUCTS

Inflammasome reporter cells

- THP1-HMGB1-Lucia™ Cells
- THP1-ASC-GFP Cells

Inflammasome inhibitors

Dectin reporter cells

- HEK-Blue™ hDectin-1a Cells
- HEK-Blue™ hDectin-1b Cells
- HEK-Blue™ mDectin-1b Cells

Reporter detection reagents

- QUANTI-Luc™ Gold

Inflammasomes : connecting innate and adaptive immunity

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
1. Martinon F. et al., 2002. The inflammasome: a molecular platform triggering activation of inflammatory caspases and processing of proIL-1 β . *Mol Cell*. 10:417-26. 2. Broz P. and Dixit V.M., 2016. Inflammasomes: mechanism of assembly, regulation and signalling. *Nat Rev Immunol*. 16:407-20. 3. Kovacs S.B. and Miao E.A., 2017. Gasdermins: effectors of pyroptosis. *Trends Cell Biol*. 27:673-84. 4. Baroja-Mazo A. et al., 2014. The NLRP3 inflammasome is released as a particulate danger signal that amplifies the inflammatory response. *Nat Immunol*. 15:738-48. 5. Broz P. et al., 2010. Differential requirement for caspase-1 autoproteolysis in pathogen-induced cell death and cytokine processing. *Cell Host Microbe*. 8:471-83. 6. Evavold C.L. et al., 2018. The pore-forming protein gasdermin D regulates interleukin-1 secretion from living macrophages. *Immunity*. 48:1-10. 7. Evavold C.L. and Kagan J.C., 2018. How inflammasomes inform adaptive immunity. *J Mol Biol*. 430:217-37. 8. Zhu J. et al., 2010. Differentiation of effector CD4 T cell populations. *Annu Rev Immunol*. 28:445-489. 9. Tang J. et al., 2018. Regulation of C-type lectin receptor-mediated antifungal immunity. *Front Immunol*. 9:123. 10. Geijtenbeek T.B.H. and Gringhuis S.I., 2016. C-type lectin receptors in the control of T helper cell differentiation. *Nat Rev Immunol*. 16:433-48. 11. Muñoz-Wolf N. et al., 2017. The role of inflammasomes in adjuvant-driven humoral and cellular immune responses. *Immunopotentiators in Modern Vaccines*. 23-42. 12. Dostert C. et al., 2013. Innate and adaptive effects of inflammasomes on T cell responses. *Curr Opin Immunol*. 25:359-65. 13. Carroll E.C. et al., 2016. The vaccine adjuvant chitosan promotes cellular immunity via DNA sensor cGAS-STING-dependent induction of type I interferons. *Immunity*. 44:597-608.

InvivoGen

ASIA
Unit 106, 1F,
15W Phase 3, Hong Kong Science Park,
Pak Shek Kok,
Hong Kong

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




Visualization of ASC speck formation by fluorescence microscopy: THP1-ASC-GFP cells were primed with 1 μ g/ml LPS-EK for 3 hours, inducing the expression of the ASC::GFP fusion protein in the cytoplasm (left panel). Cells were then incubated with 250 ng/ml complexed Poly(dA:dT) and ASC speck formation was monitored over 1 to 3 hours post-activation. In most cells, only one speck forms upon inflammasome activation (arrows in right panel). Scale bar: 50 μ m.

Pyroptosis response of THP1-HMGB1-Lucia™ cells: cells were primed with 1 μ g/ml LPS-EK for 3h and then incubated with inflammasome inducers (8 μ M Nigericin, 0.5 μ g/ml complexed Poly(dA:dT) or 100 μ g/ml CPPD). Lucia luciferase activity and LDH release in the supernatant were quantified at 2, 3, 4, 5 and 6 hours post-induction.

NF- κ B response of HEK-Blue™ IL-1 β cells: THP1-HMGB1-Lucia™ cells were incubated with inflammasome inducers (8 μ M Nigericin, 100 μ g/ml CPPD, 0.5 μ g/ml complexed Poly(dA:dT) or 1 mg/ml Zymosan). After 24h, supernatants were incubated with HEK-Blue™ IL-1 β cells overnight and SEAP activity was assessed using QUANTI-Blue™.

PRODUCT	QUANTITY	CAT. CODE
THP1-HMGB1-Lucia™ cells	3-7 x 10 ⁶ cells	thp-hmgLuc
THP1-ASC-GFP cells	3-7 x 10 ⁶ cells	thp-ascGFP
HEK-Blue™ IL-1 β cells	3-7 x 10 ⁶ cells	hkb-il1b
CPPD	5 mg	tlrl-cppd
Nigericin	10 mg	tlrl-nig
Poly(dA:dT)	200 μ g	tlrl-patn

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
Effect of inflammasome inhibitors on IL-1β release by THP1-HMGB1-Lucia™ cells: THP1-HMGB1-Lucia™ cells were primed with LPS and treated with MSU (200 µg/ml) or Poly(dA:dT) complexes (0.5 µg/ml) in the presence of various inhibitors at 0.5 µM. After overnight incubation, cell supernatants were added to HEK-Blue™ IL-1β cells for 16 hours. IL-1β levels were determined by assessing SEAP activity in the supernatant using QUANTI-Blue™.

 www.invivogen.com/inflammasome-inhibitors

Dectin-1a and Dectin-1b reporter cells

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
NF-κB responses of HEK-Blue™ hDectin-1a and -1b and HEK-Blue™ Null I-v cells to Dectin-1 agonists: 10 µg/ml Zymozan, Laminarin, WGP soluble, 100 µg/ml WGP dispersible, or 3 x 10⁶ cells/ml HKCA. 10 ng/ml TNF-α was used as a positive control. After 24h, SEAP activity was assessed in the supernatant using QUANTI-Blue™.

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1. Heinsbroek, S.E, et al, 2006. Expression of functionally different Dectin-1 isoforms by murine macrophages. *J Immunol.* 176: 5513. **2. Smith A.J., et al, 2018.** Immunoregulatory activity of the natural product Laminarin varies widely as a result of its physical properties. *J Immunol.* 200:788.

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

www.invivogen.com/reporter-detection

Procedure

Lucia luciferase activity detection with QUANTI-Luc-based detection media: THP1-Dual™ cells were stimulated for 24h with 10 µg/ml 2'3'-cGAMP. The ISG response was assessed by determining Lucia luciferase activity in the supernatant using QUANTI-Luc™, QUANTI-Luc™ Plus, or QUANTI-Luc™ Gold (QUANTI-Luc™ Plus with QLC Stabilizer). Relative light units (RLUs) were measured over 10 minutes.

PRODUCT	QUANTITY	CAT. CODE
QUANTI-Luc™ Gold	1 pouch (25 ml)	rep-qlcg1
	2 pouches (50 ml)	rep-qlcg2
	5 pouches (125 ml)	rep-qlcg5
	1 pouch (500 ml)	rep-qlcg-500
QUANTI-Luc™	2 pouches (50 ml)	rep-qlc1
	5 pouches (125 ml)	rep-qlc2
	1 pouch (500 ml)	rep-qlc-500

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www.invivogen.com/selective-antibiotics

InvivoGen

Europe Tel: +33 562 71 69 39
USA Tel: +1 888 457 5873
Asia Tel: +852 3622 3480

Fax: +33 562 71 69 30
Fax: +1 858 457 5843
Fax: +852 3622 3483

info.eu@invivogen.com
info@invivogen.com
info.hk@invivogen.com