SUMMARY :

REVIEW

RIG-I and cancer immunotherapy

PRODUCTS

Nucleic acid signaling reporter cell lines

- THP1-Dual™ KO-TBK1 Cells
- THP1-Dual™ KO-IRF3 Cells

Nucleic acid-signaling inhibitors

- BX795
- MRT67307
- RU.521
- H-151

NATE™ Nucleic acid transfection enhancer

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
1. Iurescia S. *et al.*, 2018. Nucleic acid sensing machinery: targeting innate immune system for cancer therapy. *Recent Pat. Anticancer Drug Discov.* 13:2. 2. Heidegger S. *et al.*, 2019. RIG-I activating immunostimulatory RNA boosts the efficacy of anticancer vaccines and synergizes with immune checkpoint blockade. *EBioMedicine*. 41:146. 3. Poeck H. *et al.*, 2008. 5'-triphosphate-siRNA: turning gene silencing and RIG-I activation against melanoma. *Nat. Med.* 14:1256. 4. Elion DL. *et al.*, 2018. Harnessing RIG-I and intrinsic immunity in the tumor microenvironment for therapeutic cancer treatment. *Oncotarget*. 9:29007. 5. Ren X. *et al.*, 2019. RIG-I selectively discriminates against 5'-monophosphate RNA. *Cell Reports*. 26:2019. 6. Chen ST. *et al.*, 2019. NLRP12 regulates anti-viral RIG-I activation via interaction with TRIM25. *Cell Host Microbe*. 25:602. 7. Zevini A. *et al.*, 2017. Crosstalk between cytoplasmic RIG-I and STING sensing pathways. *Trends Immunol.* 38:194. 8. Perry AK. *et al.*, 2004. Differential requirement for TANK-binding kinase-1 in type I interferon responses to Toll-like receptor activation and viral infection. *J. Exp. Med.* 199:1651. 9. Fang R. *et al.*, 2017. MAVS activates TBK1 and IKKε through TRAFs in NEMO dependent and independent manner. *PLoS Pathog.* 13(11):e1006720.

InvivoGen

ASIA

Unit 106, 1F,

15W Phase 3, Hong Kong Science Park,
Pak Shek Kok,
Hong Kong

T: +852 3622 3480
F: +852 3622 3483

info.hk@invivogen.com
www.invivogen.com

Nucleic acid signaling reporter cell lines

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Validation of KO-TBK1 and KO-IRF3 Cells


Figure 1: 2'3'-cGAMP-induced IRF or NF- κ B responses in THP-1 Dual™ derived cells. 2×10^5 THP1-Dual™ (WT), KO-TBK1 or KO-IRF3 cells were incubated with 30 μ g/ml 2'3'-cGAMP (STING agonist). After overnight incubation, ISG and NF- κ B induction were assessed by measuring the activity of Lucia luciferase and SEAP in the cell supernatant using QUANTI-Luc™ and QUANTI-Blue™ Solution, respectively. Activity fold increase over non-induced cells (Lucia luciferase readout) or reading of optical density (OD) at 630 nm (SEAP readout) are shown.

Nucleic acid signaling reporter cell collection

- RAW-Lucia ISG™ KO-cGAS
- RAW-Lucia ISG™ KO-IFI16
- RAW-Lucia ISG™ KO-IRF3
- RAW-Lucia ISG™ KO-IRF5
- RAW-Lucia ISG™ KO-IRF7
- RAW-Lucia ISG™ KO-MAVS
- RAW-Lucia ISG™ KO-MDA5
- RAW-Lucia ISG™ KO-RIG-I
- RAW-Lucia ISG™ KO-TBK1
- RAW-Lucia ISG™ KO-TREX1
- THP-1 Dual™ KO-cGAS
- THP-1 Dual™ KO-IFI16
- THP-1 Dual™ KO-MyD
- THP-1 Dual™ KO-STING
- THP-1 Dual™ KO-TREX1

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RELATED PRODUCTS

PRODUCT	DESCRIPTION	CAT. CODE
2'3'-cGAMP	mammalian cyclic dinucleotide	tlrl-nacga23-02
3p-hpRNA	5' triphosphate hairpin RNA	tlrl-hprna
5'ppp-dsRNA	5' triphosphate double stranded RNA	tlrl-3prna

www.invivogen.com/ko-cell-lines

Differential responses depending on the RIG-I agonist and transfection reagent

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
1. Perry AK. *et al.*, 2004. Differential requirement for TANK-binding kinase-1 in type I interferon responses to Toll-like receptor activation and viral infection. *J. Exp. Med.* 199:1651.
2. Deng W. *et al.*, 2008. Negative regulation of virus-triggered IFN- β signaling pathway by alternative splicing of TBK1. *J. Biol. Chem.* 283:35590.

Nucleic acid signaling inhibitors

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1. Clark K. *et al.*, 2009. Use of the pharmacological inhibitor BX795 to study the regulation and physiological roles of TBK1 and IκB kinase ε. *J. Biol. Chem.* 284:21. 2. Clark K. *et al.*, 2011. Novel cross-talk within the IKK family controls innate immunity. *Biochem. J.* 434:93. 3. Vincent J. *et al.*, 2017. Small molecule inhibition of cGAS reduces interferon expression in primary macrophages from autoimmune mice. *Nat. Commun.* 28:750. 4. Zhou W. *et al.*, 2018. Structure of the human cGAS-DNA complex reveals enhanced control of immune surveillance. *Cell.* 174:300. 5. Haag S.M. *et al.*, 2018. Targeting STING with covalent small-molecule inhibitors. *Nature.* 559:269.

BX795 and MRT67307 inhibition of TBK1/IKKε signaling


Figure 3: Inhibition of 3p-hpRNA-induced responses in THP1-Dual™ Cells: 2x10⁵ THP1-Dual™ cells were transfected with 1 μg/ml 3p-hpRNA/Lyovec™ with or without 1 μg/ml BX795 or 3 μg/ml MRT67307. After overnight incubation, activation of IRF or NF-κB was assessed by measuring Lucia luciferase or SEAP activity in the supernatant.

RU.521 and H-151 inhibition of cGAS and STING signaling


Figure 4: Inhibition of DNA-induced responses in RAW-Lucia™ ISG Cells: 10⁵ RAW-Lucia™ ISG cells were transfected with 1 μg/ml G3-YSD/LTX or incubated with 20 μg/ml 2'3'-cGAMP in the presence or absence of 2.5 μg/ml RU.521 or 130 ng/ml H-151, respectively. After overnight incubation, IRF activation was assessed by measuring Lucia luciferase activity in the supernatant.

PRODUCT	TARGET	QUANTITY	CAT. CODE
BX795	TBK1/IKKε	5mg	tlrl-bx7
H-151	STING	10 mg	inh-h151
MRT67307	TBK1/IKKε	10 mg	inh-mrt
RU.521	cGAS	2 mg	inh-ru521

 www.invivogen.com/prr-signaling-inhibitors


Figure 2: RNA induction of ISG or NF-κB in RAW-Lucia™ ISG or THP1-Dual™-derived cells. 10⁵ RAW-Lucia™ ISG, KO-TBK1 or KO-IRF3 cells (a, b) or 2x10⁵ THP1-Dual™, KO-TBK1 or KO-IRF3 cells (c-d) were transfected with 1 μg/ml 3p-hpRNA or 5'ppp-dsRNA complexed with Lyovec™ (a) or LTX (b-d). After overnight incubation, ISG induction was assessed by measuring bioluminescent activity of the Lucia luciferase in the supernatant using QUANTI-Luc™. Activity fold increase over non-transfected cells is shown (a-c). The NF-κB activity in THP1-Dual™-derived cells was assessed by measuring the SEAP activity in the supernatant using QUANTI-Blue™ Solution. Reading of optical density (OD) at 630 nm is shown (d).

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
NATE™ Nucleic Acid Transfection Enhancer NEW

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Obstacles to successful nucleic acid transfection

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
PRODUCT	QUANTITY	CAT. CODE
NATE™	1 ml (~100 reactions)	lyec-nate


www.invivogen.com/nate

Increased DNA expression in NATE™ treated cells

NATE™ enhances transient DNA expression in human THP-1 monocytes: Transfection of a ~3 kb GFP-expressing plasmid into THP-1 cells was performed using GeneXplus without (a; top; b; yellow) or with (a; bottom; b; green) the NATE™ reagent. After 48 hours, cells were visualized by fluorescence microscopy (a), and transfection efficiency was measured using flow cytometry (b). Data are presented as a fold change normalized to the transfection efficiency without NATE™ reagent.

Increased rate of stably transfected cells with NATE™

NATE™ increases the number of stable SEAP-expressing clones in murine RAW 264.7 macrophages: Transfection of a ~10 kb SEAP-expressing plasmid into RAW 264.7 cells was performed using Lipofectamine® LTX, without (left) or with (right) the NATE™ reagent. After 10 days in selection with Blasticidin, the number of stables clones expressing SEAP (blue wells) was readily visualized using QUANTI-Blue™ Solution detection reagent.

Examples of transfections with NATE™

	THP-1	RAW 264.7
Plate size	12-well	24-well
Seeding Cell Density	5 x 10 ⁵ c/w/ml	2 x 10 ⁵ c/w/ml
Transfection reagent used	GeneXplus	Lipofectamine® LTX
Volume of 100x NATE™	10 µl	10 µl

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