## SUMMARY :

### REVIEW

#### AhR's key role in the intestinal microbiota and immunity

### PRODUCTS

#### AhR Agonists

- ITE
- L-Kynurenine
- FICZ

#### AhR Inhibitor

- CH-223191

#### AhR Reporter Cell Lines

- HT29-Lucia™ AhR cells
- HepG2-Lucia™ AhR cells


#### Gut Microbiota-Related Cell Lines

#### QUANTI-Blue™ Solution

# AhR's key role in the intestinal microbiota and immunity

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
and intestinal immunity. *Mucosal Immunol.* 11(4):1024-38. 3. Gutiérrez-Vasquez C. et al., 2018. Regulation of the immune response by the aryl hydrocarbon receptor. *Immunity.* 48:19-33. 4. Hubbard T.D. et al., 2015. Indole and tryptophan metabolism: endogenous and dietary routes to Ah receptor activation. *Drug Metab. Dispos.* 43:1522-35. 5. Murray I.A. et al., 2016. Expression of the aryl hydrocarbon receptor contributes to the establishment of intestinal microbial community structure in mice. *Sci. Rep.* 6:33969. 6. Cervantes-Barragan L. & Colonna M., 2018. AHR signaling in the development and function of intestinal immune cells and beyond. *Semin. Immunopathol.* 40(4):371-77. 7. Ehrlich A.K. et al., 2018. TCDD, FICZ, and other high affinity AhR ligands dose-dependently determine the fate of CD4+ T cell differentiation. *Toxicol. Sci.* 161(2):310-20. 8. Boule L.A. et al., 2018. Aryl hydrocarbon receptor signaling modulates antiviral immune responses: ligand metabolism rather than chemical source is the stronger predictor of outcome. *Sci. Rep.* 8:1826. 9. Bessede A. et al., 2014. Aryl hydrocarbon receptor control of a disease tolerance defence pathway. *Nature.* 511:184-90. 10. Zhu J. et al., 2018. Aryl hydrocarbon receptor promotes IL-10 expression in inflammatory macrophages through Src-STAT3 signaling pathway. *Front. Immunol.* 9:2033. 11. Gao J. et al., 2018. Impact of the gut microbiota on intestinal immunity mediated by tryptophan metabolism. *Front. Cell. Infect. Microbiol.* 8:13.

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

1. Lamas B. et al., 2018. Aryl hydrocarbon receptor and intestinal immunity. *Mucosal Immunol.* 11(4):1024-38. 2. Murray I.A. et al., 2017. Ligand activation of the Ah receptor contributes to gastrointestinal homeostasis. *Curr. Opin. Toxicol.* 2:15-23. 3. Hubbard T.D. et al., 2015. Indole and tryptophan metabolism: endogenous and dietary routes to Ah receptor activation. *Drug Metab. Dispos.* 43:1522-35.

[www.invivogen.com/AhR-ligands](http://www.invivogen.com/AhR-ligands)

**Inhibitor of AhR activity by CH-223191 in HepG2-Lucia™ AhR cells:** Cells were stimulated with 30 nM TCDD or 30 μM ITE in the presence of 10 μM CH-223191. After overnight incubation, the AhR response was assessed by determining Lucia luciferase activity in the supernatant using QUANTI-Luc™. Percentages of the maximal response for the ligand with no inhibitor are shown.

1. Kim, S.H. et al., 2006. Novel compound 2-methyl-2H-pyrazole-3-carboxylic acid (2-methyl-4-o-tolylazo-phenyl)-amide (CH-223191) prevents 2,3,7,8-TCDD-induced toxicity by antagonizing the aryl hydrocarbon receptor. *Mol. Pharmacol.* 69:1871-78. 2. Zhao B. et al., 2010. CH223191 is a ligand-selective antagonist of the Ah (Dioxin) receptor. *Toxicol. Sci.* 117:393-403.

[www.invivogen.com/AhR-inhibitors](http://www.invivogen.com/AhR-inhibitors)

# Aryl Hydrocarbon Receptor Reporter Cells

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


## HT29-Lucia™ AhR Cells

**Induction of AhR activity by tryptophan byproducts in HT29-Lucia™ AhR cells:** cells were incubated with 100 μM ITE, 18 μM FICZ, or 150 μM L-Kynurenine. After overnight incubation, the AhR activation was assessed by determining Lucia luciferase activity in the supernatant using QUANTI-Luc™. Data are expressed as fold responses as compared to non-induced cells.

## HepG2-Lucia™ AhR Cells

**Induction of AhR activity by tryptophan byproducts in HepG2-Lucia™ AhR cells:** cells were incubated with 30 μM ITE, 18 μM FICZ, or 150 μM L-Kynurenine. After overnight incubation, the AhR activation was assessed by determining Lucia luciferase activity in the supernatant using QUANTI-Luc™. Data are expressed as fold responses as compared to non-induced cells.

PRODUCT	QUANTITY	CAT. CODE
HepG2-Lucia™ AhR cells	<b>NEW</b> 3-7 x 10 <sup>6</sup> cells	hpgl-ahr
HT29-Lucia™ AhR cells	<b>NEW</b> 3-7 x 10 <sup>6</sup> cells	ht2l-ahr
FICZ	<b>NEW</b> 1 mg	tlrl-ficz
ITE	<b>NEW</b> 10 mg	tlrl-ite
L-Kynurenine	<b>NEW</b> 10 mg	tlrl-kyn
CH-223191	<b>NEW</b> 10 mg	inh-ch22
Quanti-Luc™	2 pouches (2 x 25 ml)	rep-qlc1
Quanti-Luc™ Gold	1 pouch (25 ml)	rep-qlcg1
Zeocin™	5 x 1 ml	ant-zn-05

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- TLR Reporter Cells
- NOD1/2 Reporter Cells
- Dectin-1a/b Reporter Cells

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Please see our website for full lists of human and murine PRR reporter cell lines and their agonists.

1. Swiatczak, B. & Cohen, I.R. 2015. Gut feelings of safety: tolerance to the microbiota mediated by innate immune receptors. *Microbiol Immunol.* 59:573-585. 2. Mu, C. et al. 2015. Crosstalk between the immune receptors and gut microbiota. *Curr Protein Pept Sci.* 16:622-631.

 [www.invivogen.com/reporter-cells](http://www.invivogen.com/reporter-cells)

## QUANTI-Blue™ Solution NEW

### Liquid formulation for SEAP detection

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
PRODUCT		QUANTITY	CAT. CODE
QUANTI-Blue™ Solution	<b>NEW</b>	5 ml	rep-qbs
QUANTI-Blue™ Solution	<b>NEW</b>	Bulk	Please enquire

 [www.invivogen.com/quant-blue](http://www.invivogen.com/quant-blue)

**InvivoGen**

Europe Tel: +33 562 71 69 39  
USA Tel: +1 888 457 5873  
Asia Tel: +852 3622 3480

Fax: +33 562 71 69 30  
Fax: +1 858 457 5843  
Fax: +852 3622 3483

info.eu@invivogen.com  
info@invivogen.com  
info.hk@invivogen.com