**SUMMARY :**

**REVIEW**

**Interferons λs: guardians of the front-lines**

**PRODUCTS**

**Interferon-λ Reporter Cell Line**

- HEK-Blue™ IFN-λ cells

**Interferon-λ Antibodies**

- Anti-hIL-28a-IgG (hIFN-λ2)
- Anti-hIL-28b-IgG (hIFN-λ3)
- Anti-hIL-29-IgG (hIFN-λ1)

**JAK/STAT Signaling Inhibitors**

- CP-690550
- CYT387
- Ruxolitinib

**Recombinant Type I Interferons**

- Recombinant human interferon-α

**Mycoplasma Detection & Elimination**

# Interferon λs: guardians of the front-lines

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
1. Pestka S., 2007. The Interferons: 50 years after their discovery there is much more to learn. *J. Biol. Chem.* 282: 20047. 2. Kotenko S.V. and Durbin J.E., 2017. Contribution of type III interferons to antiviral immunity: location, location, location. *J. Biol. Chem.* 292: 7295. 3. Lazear H.M. et al., 2015. Interferon-λ: Immune Functions at Barrier Surfaces and Beyond. *Immunity* 43: 15. 4. Durbin R.K. et al., 2013. Diverse intracellular pathogens activate Type III interferon expression from peroxisomes. *Nat. Immunol.* 15:717. 5. Odendall C. et al., 2014. Dynamic expression profiling of type I and type III interferon-stimulated hepatocytes reveals a stable hierarchy of gene expression. *Hepatology*, 59: 1262. 6. Andreakos E. et al., 2017. Interferon-λs: front-line guardians of immunity and homeostasis in the respiratory tract. *Front. Immunol.* 8: 1232. 7. Syedbasha M. and Egli A., 2017. Interferon Lambda: Modulating Immunity in Infectious Diseases. *Frontiers in Immunology* 8: 119. 8. Obajemu A.A. et al., 2017. IFN-λ4 Attenuates Antiviral Responses by Enhancing Negative Regulation of IFN Signaling. *J. Immunol.* 199: 3808. 9. Zanon I. et al., 2017. Interferon (IFN)-λ Takes the Helm: Immunomodulatory Roles of Type III IFNs. *Front. Immunol.* 8: 1661.

# Interferon Lambda Reporter Cell Line

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
**Response of HEK-Blue™ IFN-λ cells to type I, II and III IFNs:** HEK-Blue™ IFN-λ cells were incubated with increasing concentrations of recombinant human or mouse IFN-λ (ng/ml) or human IFN-α, IFN-β or IFN-γ (IU/ml). After 24h incubation, ISG activation was assessed by measuring SEAP levels in the supernatant using Quanti-Blue™. EC<sub>50</sub> is indicated for each cytokine (N/A: non applicable).


[www.invivogen.com/cytokine-reporter-cells](http://www.invivogen.com/cytokine-reporter-cells)

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**Neutralizing activity of interferon lambda antibodies:** HEK-Blue™ IFN-α/β cells were incubated with 10 ng/ml recombinant hIL-28a, hIL28b, hIL-29 and increasing concentrations of their cognate antibody. After 24h incubation, recombinant IFN-λ induced ISG activation was assessed by measuring SEAP levels in the supernatant using Quanti-Blue™. Percentages of maximal response (no antibody) and IC<sub>50</sub> for 10 ng/ml of each cytokine is shown.

[www.invivogen.com/antibodies](http://www.invivogen.com/antibodies)

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


1. Roskoki R, Jr, 2016. Janus kinase (JAK) inhibitors in the treatment of inflammatory and neoplastic diseases. *Pharmacol. Res.* 111:784-803.

 [www.invivogen.com/jakstat-inhibitors](http://www.invivogen.com/jakstat-inhibitors)

Effect of JAK/STAT inhibitors on HEK-Blue™ IFN- $\alpha/\beta$  cell response to type I and type III IFNs: HEK-Blue™ IFN- $\alpha/\beta$  cells were incubated with 3 U/ml hIFN- $\alpha$ 2b (grey), 1 U/ml hIFN- $\beta$ 1 (purple) or 10 ng/ml hIL-29 (hIFN- $\lambda$ 1) (red) and increasing concentrations of JAK inhibitors. After 24h incubation, IFN-induced ISG activation was assessed by measuring SEAP levels in the supernatant using Quanti-Blue™. Percentages of maximal response (no inhibitor) for each cytokine are shown.

PRODUCT	QUANTITY	CAT. CODE
CP-690550 (Tofacitinib)	5 mg	tlrl-cp69
CYT387 (Mometinib)	10 mg	inh-cy87
Ruxolitinib (INC424)	5 mg	tlrl-rux

# Recombinant Type I Interferons


## Recombinant human interferon- $\alpha$ s NEW

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1. Hoffmann H-H, et al., 2015. Interferons and viruses: an evolutionary arms race of molecular interactions. *Trends Immunol.* 36:124. 2. Schreiber G. & Piehler J., 2015. The molecular basis for functional plasticity in type I interferon signaling. *Trends Immunol.* 36:139. 3. Kurunganti S, et al., 2014. Production and characterization of thirteen human type-I interferon- $\alpha$  subtypes. *Protein Expr. Purif.* 103:75.

Response of THP1-Dual™ cells to recombinant human IFN- $\alpha$ s: THP1-Dual™ cells were incubated with increasing concentrations of each recombinant human IFN- $\alpha$  subtype (ng/ml). After 24h incubation, ISG activation was assessed by measuring Lucia luciferase activity in the supernatant using Quanti-Luc™.

PRODUCT	QUANTITY	CAT. CODE
Recombinant hIFN- $\alpha$ "n"	1 $\mu$ g	rcyc-hifna"n"
Recombinant hIFN- $\alpha$ set	12 x 1 $\mu$ g	rcyc-hifna
THP1-Dual™ cells	3-7 x 10 <sup>6</sup> cells	thpd-nfs

"n" refers to the number of the IFN- $\alpha$  subtype, e.g. rhIFN- $\alpha$ 2.

Recombinant human IFN- $\alpha$ s comprise IFN- $\alpha$ 1 (D), - $\alpha$ 2 (2b), - $\alpha$ 4 (4a, M1), - $\alpha$ 5 (G), - $\alpha$ 6 (K), - $\alpha$ 7 (J1), - $\alpha$ 8 (B2), - $\alpha$ 10 (C), - $\alpha$ 14 (H2), - $\alpha$ 16 (WA), - $\alpha$ 17 (I) and - $\alpha$ 21 (F).

 [www.invivogen.com/human-ifna](http://www.invivogen.com/human-ifna)

# Mycoplasma Detection and Elimination

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1. Drexler H.G and Uphoff C.C., 2002. Mycoplasma contamination of cell cultures: incidence, sources, effects, detection, elimination, prevention. *Cytotechnology*, 39:75. 2. Rottem S. 2003. Interaction of Mycoplasmas with host cells. *Physiol Rev*, 83:417. 3. Zakharova E. et al., 2010. Mycoplasma suppression of THP-1 Cell TLR responses is corrected with antibiotics. *PLoS One*, 5(3):e9900.

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

- ### TOP 5 REASONS TO TEST
1. 20 to 35% worldwide contamination, including cell banks
  2. Serious impact on data reliability and reproducibility
  3. Loss of cell lines
  4. Financial impact
  5. Mycoplasma testing required by most journals

PRODUCT	QUANTITY	CAT. CODE
PlasmoTest™	1 kit (250 tests)	rep-pt1
PlasmoTest™ controls	200 tests	pt-ctr2
PlasmoTest™ refills	500 tests	rep-ptrk
Plasmocin™	25 mg (1 ml)	ant-mpt-1
Plasmocure™	100 mg (1 ml)	ant-pc

[www.invivogen.com/mycoplasma](http://www.invivogen.com/mycoplasma)