WINTER 2017/2018

SUMMARY :

REVIEW

The RIG-I and STING Alliance

PRODUCTS

RIG-I Pathway

- RIG-I agonists
- RIG-I/MAVS reporter cells

STING Variants

- SAVI STING reporter cells
- Recombinant human STING

Nucleic Acid Complexing Agent

- LyoVec™

Selective Antibiotics

- Cell culture-grade antibiotics

InvivoGen

ASIA

Unit 106, 1F,
15W Phase 3, Hong Kong Science Park,
Pak Shek Kok,
Hong Kong

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
1. Gebhardt A. et al., 2017. Discrimination of Self and Non-Self Ribonucleic Acids. *Journal of Interferon & Cytokine Research* 37: 184-97.
2. Tao J. et al., 2016. cGAS-cGAMP-STING: The three musketeers of cytosolic DNA sensing and signaling. *IUBMB Life* 68: 858-70.
3. Cheng G. et al., 2007. Double-stranded DNA and double-stranded RNA induce a common antiviral signaling pathway in human cells. *Proceedings of the National Academy of Sciences* 104: 9035-40.
4. Ablasser A. et al., 2009. RIG-I-dependent sensing of poly(dA:dT) through the induction of an RNA polymerase III-transcribed RNA intermediate. *Nat Immunol* 10: 1065-72.
5. Chiu Y.-H. et al., 2009. RNA Polymerase III Detects Cytosolic DNA and Induces Type I Interferons through the RIG-I Pathway. *Cell* 138: 576-91.
6. Choi M.K. et al., 2009. A selective contribution of the RIG-I-like receptor pathway to type I interferon responses activated by cytosolic DNA. *Proceedings of the National Academy of Sciences* 106: 17870-5.
7. Zhong B. et al., 2008. The Adaptor Protein MITA Links Virus-Sensing Receptors to IRF3 Transcription Factor Activation. *Immunity* 29: 538-50.
8. Ishikawa H. & Barber GN. 2008. STING is an endoplasmic reticulum adaptor that facilitates innate immune signalling. *Nature*. 455(7213):674-8.
9. Castanier C. et al., 2010. Mitochondrial dynamics regulate the RIG-I-like receptor antiviral pathway. *EMBO reports* 11: 133-8.
10. Gao D. et al., 2013. Cyclic GMP-AMP Synthase Is an Innate Immune Sensor of HIV and Other Retroviruses. *Science* 341: 903-6.
11. Mankan AK. et al., 2014. Cytosolic RNA:DNA hybrids activate the cGAS-STING axis. *EMBO J.* 33(24):2937-46.
12. Civril F. et al., 2013. Structural mechanism of cytosolic DNA sensing by cGAS. *Nature* 498: 332-7.
13. Liu Y. et al., 2016. RIG-I-Mediated STING Upregulation Restricts Herpes Simplex Virus 1 Infection. *Journal of Virology* 90: 9406-19.
14. Wu X. et al., 2017. RIG-I and IL-6 are negative-feedback regulators of STING induced by double-stranded DNA. *PLOS ONE* 12: e0182961.

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www.invivogen.com/rlr-cell-lines

ISG Responses to RIG-I ligands: A549-derived cells were stimulated with 300 ng/ml 5'ppp-dsRNA or 30 ng/ml 3p-hpRNA, while RAW- and HEK-derived cells were induced with 1 µg/ml 5'ppp-dsRNA or 3p-hpRNA. After overnight incubation, the ISG response was assessed by determining Lucia luciferase activity in the supernatant using QUANTI-Luc™ expressed as relative light units (RLUs).

PRODUCT	QUANTITY	CAT. CODE
A549-Dual™ cells	3-7 x 10 ⁶ cells	a549d-nfis
A549-Dual™ KO-MAVS cells	3-7 x 10 ⁶ cells	a549d-komavs
A549-Dual™ KO-RIG-I cells	3-7 x 10 ⁶ cells	a549d-korigi
RAW-Lucia™ ISG cells	3-7 x 10 ⁶ cells	rawl-isg
RAW-Lucia™ ISG-KO-MAVS cells	3-7 x 10 ⁶ cells	rawl-komavs
RAW-Lucia™ ISG-KO-RIG-I cells	3-7 x 10 ⁶ cells	rawl-korigi
HEK-Lucia™ Null cells	NEW 3-7 x 10 ⁶ cells	hkl-null
HEK-Lucia™ RIG-I cells	NEW 3-7 x 10 ⁶ cells	hkl-hrigi
5'ppp-dsRNA	25 µg	tlrl-3prna
3p-hpRNA	NEW 25 µg	tlrl-hprna

RELATED PRODUCTS

PRODUCT	DESCRIPTION	CAT. CODE
QUANTI-Luc™	Luciferase detection reagent	rep-qlc1
QUANTI-Blue™	SEAP detection reagent	rep-qb1

For more information, please visit : www.invivogen.com

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


1. Munoz J. et al., 2015. Stimulator of interferon genes-associated vasculopathy with onset in infancy: a mimic of childhood granulomatosis with polyangiitis. *JAMA Dermatology* 151: 872-7.
2. Jeremiah N. et al., 2014. Inherited STING-activating mutation underlies a familial inflammatory syndrome with lupus-like manifestations. *The Journal of Clinical Investigation* 124: 5516-20.
3. Yi G. et al., 2013. Single nucleotide polymorphisms of human STING can affect innate immune response to cyclic dinucleotides. *PLOS ONE* 8: e77846.

PRODUCT		QUANTITY	CAT. CODE
THP1-Dual™ KI-hSTING-S154 (SAVI) cells	NEW	3-7 x 10 ⁶ cells	thpd-s154
THP1-Dual™ KI-hSTING-M155 (SAVI) cells	NEW	3-7 x 10 ⁶ cells	thpd-m155

RELATED PRODUCTS

PRODUCT	DESCRIPTION	CAT. CODE
THP1-Dual™ KI-hSTING-R232 cells	R232 variant- ("wild-type") expressing cells	thpd-r232
2'3'-cGAMP	Cyclic dinucleotide	tlrl-nacga23
2'3'-c-di-AM(PS)2 (Rp,Rp)	Cyclic dinucleotide	tlrl-nacda2r-01


ISG Responses to IFN-β or STING ligands: THP1-Dual™-derived cells were stimulated with 10⁴ U/ml IFN-β, 10 μg/ml 2'3'-cGAMP or 10 μg/ml 2'3'-cGAMP(PS)₂(Rp/Sp), 30 μM BX795 or 10 μg/ml Brefeldin A. After overnight incubation, the ISG response was assessed by determining Lucia luciferase activity in the supernatant using QUANTI-Luc™. Bars represent the % activity normalized on the background signal in non-induced (NI) cells.

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PRODUCT		QUANTITY	CAT. CODE
Recombinant Human STING	NEW	25 μg	rec-hsting

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
PRODUCT	QUANTITY	CAT. CODE
LyoVec™	10 ml (200 reactions)	lyec-1
LyoVec™	20 ml (400 reactions)	lyec-2

 www.invivogen.com/transfection-reagents

Procedure

01
Reconstitute lyophilized LyoVec™ with 2 ml of sterile H₂O or PBS

02
Prepare complex at 10 µg/ml:
- In a 1.5 ml tube, add first 1 µl ligand @ 1 mg/ml
- Add drop by drop 100 µl of LyoVec™
- Homogenize gently by tapping the tube
- Incubate for 15-30 minutes at room temperature

Prepare a dilution range

03
Dispense 20 µl of complex in a well of a 96-well plate containing 180 µl of cell suspension

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Selective antibiotic	Resistance gene	Concentration in bacteria	Concentration in mammalian cells
Blasticidin	<i>bsr</i> gene	25-100 µg/ml	1-10 µg/ml
G418 (Geneticin)	<i>neo</i> gene	-	400-1000 µg/ml
Hygromycin B Gold	<i>hph</i> gene	50-100 µg/ml	50-200 µg/ml
Puromycin	<i>pac</i> gene	-	1-10 µg/ml
Zeocin™	<i>Sh ble</i> gene	25 µg/ml	50-400 µg/ml

PRODUCT	QUANTITY	CAT. CODE
Blasticidin	100 mg (10 x 1 ml)	ant-bl-1
G418 (Geneticin)	1 g (10 x 1 ml)	ant-gn-1
Hygromycin B Gold	1 g (10 x 1 ml)	ant-hg-1
Puromycin	100 mg (10 x 1 ml)	ant-pr-1
Zeocin™	1 g (10 x 1 ml)	ant-zn-1

 www.invivogen.com/selective-antibiotics