CONTENTS:

EDITORIAL

The race against SARS-CoV-2 variants

PRODUCTS

Anti-Spike detection assays

- ELISA and LIPS
 - RBD-Lucia proteins
- Flow cytometry (FACS)
 - pUNO1-Spike-dfur plasmids
 - Anti-CoV2RBD-cas-mIgG2a
 - Anti-CoV2RBD-imd-mIgG2a
 - Anti-CoV2RBD-bam-mIgG2a
 - Anti-CoV2RBD-ete-mIgG2a

Neutralization Assays

- Spike-pseudotyped lentiviral infection
 - pLV-Spike plasmids
 - A549-hACE2-TMPRSS2 Cells
 - HEK-Blue™ hACE2-TMPRSS2 Cells
- Spike-ACE2 dependent cell fusion
 - pUNO1-Spike plasmids
 - 293-hMyD88 Cells
 - HEK-Blue™ hACE2 Cells
 - A549-Dual™ hACE2-TMPRSS2 Cells

Mycoplasma Detection

- MycoStrip™

InvivoGen

ASIA

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'x' refers to the in-house number of the Spike variant (see table on the front page).
More info at: www.invivogen.com/sars2-spike-vectors


Fig 1: Binding affinity of Anti-Spike antibodies. (A) ELISA was performed by first adding an Anti-mIgG antibody, then Anti-CoV2RBD-cas-mIgG2a (Casirivimab-derived mAb), and lastly RBD-Lucia, RBD-LuciaV2, RBD-LuciaV3, RBD-LuciaV4 or RBD-LuciaV7 protein. (B) LIPS was performed by mixing Protein A beads, a vaccinee serum, and a RBD-Lucia protein (same as (A)). (A+B) Binding affinity was assessed by measuring Lucia activity using QUANTI-Luc™.

PRODUCTS	QTY	CAT. CODE
RBD-LuciaV'X'	20 µg	rbd-lucia-v'x'

'x' refers to the in-house number of the Spike variant (see table on the front page).
More info at: www.invivogen.com/sarscov2-lucia

1. Ouyang, W. et al. 2021. *Int. J. Mol. Sci.* 22(14), 7540. doi: 10.3390/ijms22147540.


Fig 2: Detection of Spike variants by FACS. 293 parental (grey), 293-SARS2-S-dfur (original; pink), and 293-SARS2-S-dfur-V5 (P.1; orange) cells were incubated with either (A) Anti-CoV2RBD-ete-mIgG2a or (B) Anti-CoV2RBD-bam-mIgG2a for 1h at 4°C. Cells were then washed and incubated with a PE-conjugated goat anti-murine IgG for 1h at 4°C. Cell surface staining was analyzed by FACS.

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Spike-pseudotyped lentiviral infection

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'x' refers to the in-house number of the Spike variant (see table on the front page).

More info at: www.invivogen.com/293-hmyd88


Fig 3: Spike-pseudotyped particle-based inhibition assay. A549-hACE2-TMPRSS2 cells were cultured in the presence of Anti-CoV2RBD-cas-mlgG2a, Anti-CoV2RBD-imd-mlgG2a, Anti-CoV2RBD-bam-mlgG2a or Anti-CoV2RBD-ete-mlgG2a antibodies. Either (A) B.1.1.7 (pLV-SpikeV2) or (B) B.1.351 (pLV-SpikeV3) Spike-pseudotyped GFP lentiviral particles were added to the cells. After 72h, the transduction efficiency (% infectivity) was evaluated by FACS.

PRODUCTS	QTY	CAT. CODE
pLV-SpikeV'X'	20 µg	plv-spike-v'x'
A549-hACE2-TMPRSS2 Cells	3 - 7 x 10 ⁶ cells	a549-hace2tpsa
HEK-Blue™ hACE2-TMPRSS2 Cells	3 - 7 x 10 ⁶ cells	hkb-hace2tpsa

'x' refers to the in-house number of the Spike variant (see table on the front page).

More info at: www.invivogen.com/spike-pseudotyping-vectors


Fig 4: Cell fusion inhibition assay. 293-hMyD88 cells were transfected with either the pUNO1-Spike or pUNO1-SpikeV5 plasmid and then incubated with Anti-CoV2RBD-cas-mlgG2a, Anti-CoV2RBD-imd-mlgG2a, or Anti-βGal antibodies. These cells were then co-cultured with HEK-Blue™ hACE2 cells. After overnight incubation, cell fusion was assessed by measuring the activity of SEAP in the supernatant using QUANTI-Blue™ Solution.

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
Mycoplasma Detection with MycoStrip™

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
PRODUCT	DESCRIPTION	QTY	CAT. CODE
MycoStrip™	Mycoplasma contamination detection kit (strips)	10 tests	rep-mys-10
		20 tests	rep-mys-20
		50 tests	rep-mys-50


www.invivogen.com/mycostrip

Detection with MycoStrip™


DETECT MYCOPLASMA IN 4 EASY STEPS


1 PREPARE YOUR SAMPLE (~5 MINS)

2 INCUBATE AT 65°C FOR 40 MINUTES

3 STOP THE REACTION AND LOAD THE CASSETTE

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For more information visit: www.invivogen.com/mycoplasma-elimination

InvivoGen

Europe Tel: +33 562 71 69 39
USA Tel: +1 888 457 5873
Asia Tel: +852 3622 3480

Fax: +33 562 71 69 30
Fax: +1 858 457 5843
Fax: +852 3622 3483

info.eu@invivogen.com
info@invivogen.com
info.hk@invivogen.com