# CELL CULTURE CONTAMINATION

## OVERVIEW OF MICROBIAL CONTAMINANTS

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- Protect your cells
- Mycoplasma contaminations
- Bacterial contaminations
- Endotoxin contaminations
- Fungal contaminations

# InvivoGen *infocus*

## A PRACTICAL GUIDE TO AVOID CELL CULTURE CONTAMINATION

4  
PAGE

### **YOUR CELLS ARE PRECIOUS, PROTECT THEM!**

Detection  
Prevention  
Elimination

5  
PAGE

### **MYCOPLASMA CONTAMINATIONS**

Mycoplasma features & dangers  
Detection of mycoplasma contamination  
Elimination of mycoplasma

10  
PAGE

### **BACTERIAL CONTAMINATIONS**

The usual suspects  
How to detect bacterial contamination?  
Elimination of bacteria

12  
PAGE

### **ENDOTOXIN CONTAMINATIONS**

What are endotoxin features?  
What are the risks for my experiments?  
Detection of endotoxins in biological reagents

14  
PAGE

### **FUNGAL CONTAMINATIONS**

How to detect fungal contamination in cell cultures?  
Elimination of fungi

15  
PAGE

### **SUMMARY TABLE OF INVIVOGEN'S ANTIMICROBIAL AGENTS**

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
**NF-κB**: nuclear factor kappa-light-chain enhancer of activated B cells  
**PCR**: polymerase chain reaction  
**PRR**: pattern recognition receptor  
**rRNA**: ribosomal ribonucleic acid  
**SEAP**: secreted embryonic alkaline phosphatase  
**TLRs**: toll-like receptors

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

section. Editors reserve the right to demand that the data be removed from the paper if the justification is deemed unsatisfactory. In addition, authors must identify the source of cell lines (with catalog number) if they have been authenticated. InvivoGen authentication testing will be provided upon request. As of May 2015, there has been an increase in the issue of cell line misidentification. InvivoGen is strongly encouraged to comply with these reporting criteria. It is good practice to obtain cell lines from reputable repositories, to routinely authenticate cell line stocks and test them for mycoplasma contamination. Resources on cell line authentication follow.

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# YOUR CELLS ARE PRECIOUS, PROTECT THEM!

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
### Major impacts of mycoplasma contamination on cell functions

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# DETECTION OF MYCOPLASMA CONTAMINATION

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
Isothermal PCR

Immunochromatographic strip

PRODUCT	DESCRIPTION	QTY	CAT. CODE
MycoStrip™	Mycoplasma contamination detection kit (strips)	10 tests	rep-mys-10
		20 tests	rep-mys-20
		50 tests	rep-mys-50

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

### MycoStrip™ procedure


01 Combine: Reaction mix + Reaction buffer + Sample

02 Perform isothermal PCR

03 Add Migration buffer and transfer sample onto the cassette

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


Colorimetric read-out

### PlasmoTest™ procedure


01 Inactivate endogenous alkaline phosphatase in the sample


02 Add heated sample to HEK-Blue™-2 sensor cells cultured in HEK-Blue™ detection medium

03 Mycoplasma lipopeptide detection by HEK-Blue™-2 cells

04 SEAP reporter activity  
• Naked-eye: Purple/blue: Positive  
Pink: Negative  
• Measure OD 620-655nm (optional)

### MYCOPLASMA DETECTION METHODS

	MycoStrip™	PlasmoTest™
Target	16S rRNA gene	Lipoproteins
Method	Isothermal PCR	Reporter bioassay
Ease of use	++++	++
Specificity	At least the 8 most common species	
Sensitivity	$10 \cdot 10^2$ CFU/ml	$10^2 \cdot 10^5$ CFU/ml
Experiment duration	<1 hr	<1 hr (hands on) OVN (incubation)
Additional equipments/reagents	Heating bath/block	Incubator Spectrophotometer
Naked-eye data visualization	Yes	Yes
Cost	Reasonable	Cost-effective

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PRODUCT	DESCRIPTION	QTY	CAT. CODE
PlasmoTest™	Mycoplasma contamination detection kit (cells)	1 kit (250 tests)	rep-pt1
PlasmoTest™ Controls	Controls for PlasmoTest™ detection kit	200 tests	pt-ctr2
PlasmoTest™ Refills	Reagents for PlasmoTest™ detection kit	500 samples	rep-ptrk

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
## Tips for successful mycoplasma elimination

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
## Mycoplasma surveillance and elimination

Upon reception of a new cell line or maintenance of cell cultures

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## MYCOPLASMA ELIMINATION REAGENTS

	Plasmocin™	Plasmocure™
Treatment duration	2 weeks	2 weeks
Ease of use	++++	++++
Efficacy	+++	++++
Cytotoxicity	+/-	+/-
Resistance <sup>11</sup>	-	--

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## They trust InvivoGen


- Kapralov A. A. et al., 2020. *Nat Chem Biol*. Redox lipid reprogramming commands susceptibility of macrophages and microglia to ferroptotic death. DOI: 10.1038/s41589-019-0462-8.
- Kazemiha V. M. et al., 2019. *Cell J*. Effectiveness of Plasmocure™ to eliminate mycoplasma species from contaminated cell cultures: A comparative study versus other antibiotics. DOI: 10.22074/cellj.2019.5996

PRODUCT	DESCRIPTION	QTY	CAT. CODE
Plasmocin™ prophylactic	Reagent for preventing mycoplasma contamination	25 mg (10 x 1 ml)	ant-mpp
Plasmocin™ treatment	Mycoplasma elimination reagent	25 mg (1 x 1 ml) 50 mg (2 x 1 ml)	ant-mpt-1 ant-mpt
Plasmocure™	Mycoplasma elimination reagent	100 mg (1 ml)	ant-pc

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## Detection of bacterial contamination

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# ELIMINATION OF BACTERIA

## Contamination preventive reagent

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## Anti-microbial agent for primary cells

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PRODUCT	DESCRIPTION	QTY	CAT. CODE
Normocin™	Reagent for preventing microbial contamination	500 mg (10 x 1 ml) 1 g (1 x 20 ml)	ant-nr-1 ant-nr-2
Normocure™	Microbial contamination elimination reagent	100 mg (2 x 1 ml)	ant-noc
Primocin™	Reagent for preventing microbial contamination in primary cells	500 mg (10 x 1 ml) 1 g (1 x 20 ml)	ant-pm-1 ant-pm-2

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

## Endotoxin features


## Sources of contamination

## Risks for *in vitro* and *in vivo* experiments

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# DETECTION OF ENDOTOXINS IN BIOLOGICAL REAGENTS

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
### Principle of LPS Detection Kit 2

#### They trust InvivoGen

- Gray M. A. *et al.*, 2020. *Nat Chem Biol*. Targeted glycan degradation potentiates the anticancer immune response in vivo. DOI: 10.1038/s41589-020-0622-x.
- Imbert P. R. C. *et al.*, 2021. *Curr Biol*. An Acquired and Endogenous Glycocalyx Forms a Bidirectional "Don't Eat" and "Don't Eat Me" Barrier to Phagocytosis. DOI: 10.1016/j.cub.2020.09.082.

PRODUCT	DESCRIPTION	QTY	CAT. CODE
HEK-Blue™ LPS Detection Kit 2	Assay for the detection and quantification of biologically active LPS	1 kit	rep-lps2
HEK-Blue™ Selection	Antibiotics for maintenance of HEK-Blue™ cells	10 x 1 ml	hb-sel
HEK-Blue™ Endotoxin Standard	Standardized <i>E. coli</i> 055:B5 LPS	10 x 50 EU	rep-hbes-10
QUANTI-Blue™ Solution	Alkaline phosphatase detection medium (liquid form, 100X)	5 ml 10 ml	rep-qbs rep-qbs2

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
## How to detect fungal contamination in cell cultures?

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### They trust InvivoGen

- Busslinger G.A. et al., 2021. *Cell Rep.* Human gastrointestinal epithelia of the esophagus, stomach, and duodenum resolved at single-cell resolution. DOI: 10.1016/j.celrep.2021.108819
- Lebreton F. et al., 2019. *Nat Commun.* Insulin-producing organoids engineered from islet and amniotic epithelial cells to treat diabetes. DOI: 10.1038/s41467-019-12472-3.

## ELIMINATION OF FUNGI

### For preventive and removal treatments

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PRODUCT	DESCRIPTION	QTY	CAT. CODE
Fungin™	Reagent for preventing fungal contamination	75 mg (5 x 1.5 ml) 200 mg (1 x 20 ml)	ant-fn-1 ant-fn-2

# SUMMARY TABLE OF INVIVOGEN'S ANTIMICROBIAL AGENTS

CONTAMINANT	PRODUCT	DESCRIPTION	PAGE
Mycoplasma	MycoStrip™	Mycoplasma contamination detection kit	6
	PlasmoTest™	Mycoplasma contamination detection kit	7
	PlasmoTest™ Controls	Controls for PlasmoTest™ detection kit	7
	PlasmoTest™ Refills	Reagents for PlasmoTest™ detection kit	7
	Plasmocin™ prophylactic	Reagent for preventing mycoplasma contamination	8
	Plasmocin™ treatment	Mycoplasma elimination reagent	8
	Plasmocure™	Mycoplasma elimination reagent	8
Bacteria	Normocin™	Reagent for preventing microbial contamination	11
	Normocure™	Microbial contamination elimination reagent	11
	Primocin™	Reagent for preventing microbial contamination in primary cells	11
Endotoxins	HEK-Blue™ LPS Detection Kit 2	Assay for the detection and quantification of biologically active LPS	13
	HEK-Blue™ Selection	Antibiotics for maintenance of HEK-Blue™ cells	13
	HEK-Blue™ Endotoxin Standard	Standardized <i>E. coli</i> 055:B5 LPS	13
	QUANTI-Blue™ Solution	Alkaline phosphatase detection medium (liquid form)	13
Fungi	Fungin™	Reagent for preventing and eliminating fungal contamination	14

## REFERENCES

- Drexler H.G and Uphoff C.C, 2002. Mycoplasma contamination of cell cultures: incidence, sources, effects, detection, elimination, prevention. *Cytotechnology*. 39:75.
- Armstrong E.E. *et al.*, 2010. The scope of mycoplasma contamination within the biopharmaceutical industry. *Biologicals*. 38(2):211-3.
- Takeuchi O. *et al.*, 2000. Cutting Edge: Preferentially the R-Stereoisomer of the Mycoplasma Lipopeptide Macrophage-Activating Lipopeptide-2 Activates Immune Cells Through a Toll-Like Receptor 2- and MyD88- Dependent Signaling Pathway. *J. Immunol.* 164:554-557.
- Miller C.J. *et al.*, 2003. Mycoplasma infection significantly alters microarray gene expression profiles. *Biotechniques*. 35(4):812-4
- Rottem S. 2003. Interaction of Mycoplasmas with host cells. *Physiol Rev.* 83:417.
- Zakharova E. *et al.*, 2010. Mycoplasma suppression of THP-1 Cell TLR responses is corrected with antibiotics. *PLoS One*. 5(3):e9900
- Molla Kazemiha V *et al.*, 2011. Efficiency of Plasmocin™ on various mammalian cell lines infected by mollicutes in comparison with commonly used antibiotics in cell culture: a local experience. *Cytotechnology*. Dec;63(6):609-20
- Uphoff CC *et al.*, 2012. Treatment of mycoplasma contamination in cell cultures with Plasmocin. *J Biomed Biotechnol.* 2012:26767.
- Rongvaux A *et al.*, 2014. Development and function of human innate immune cells in a humanized mouse model. *Nat Biotechnol.* 32(4):364-72.
- Baronti C. *et al.*, 2013. Mycoplasma removal: simple curative methods for viral supernatants. *J Virol Methods.* 187(2):234-7.
- Deng F. *et al.*, 2012. Generation of induced pluripotent stem cells from human Tenon's capsule fibroblasts. *Mol Vis.* 18:2871-81.
- Molla Kazemiha V. *et al.*, 2019. Effectiveness of Plasmocure in elimination of mycoplasma species from contaminated cell cultures: a comparative study versus other antibiotics. *Cell. J.* 21(2):143-149.
- Romorini L *et al.*, 2013. Effect of antibiotics against Mycoplasma sp. on human embryonic stem cells undifferentiated status, pluripotency, cell viability and growth. *Plos One.* Jul 30;8(7):e70267
- Ryan J., 2008. Understanding and managing cell culture contamination. Corning Life Sciences, technical literature.
- Lincoln C. & Gabridge M. 1998. Cell culture contamination: sources, consequences, prevention, and elimination. *Methods in cell biology.* 57:49-65.
- Fogh J., 1973. Contamination in Tissue Culture, published by Academic Press Inc. 4. McGarrity G.J. & Coriell LL., 1971. Procedures to reduce contamination of cell cultures. *In Vitro.* 6(4):257-65.
- McGarrity, G. J. *et al.*, 1984. Cytogenetic effects of mycoplasma infection of cell cultures: a review. *In vitro.* 20(1):1-18.
- Gray JS. *et al.*, 2010. Got black swimming dots in your cell culture? Identification of *Achromobacter* as a novel cell culture contaminant. *Biologicals.* 38(2):273-7
- McGowan JE Jr., 2006. Resistance in nonfermenting gram-negative bacteria: multidrug resistance to the maximum. *Am J Med.* 119(6 Suppl 1):S29-36; discussion S62-70.
- Gorbet MB. & Sefton MV., 2005. Endotoxin : the uninvited guest. *Biomaterials.* 26:6811-6817.
- Leider, R. *et al.*, 2013. Endotoxin - The invisible companion in biomaterial research. *Tissue Eng. Part B. Rev.* 19(5):391-402.
- Nims, RW. & Price P.J. 2017. Best practices for detecting and mitigating the risk of cell culture contaminants. *In vitro Cell. Dev. Biol. Anim.* 53(10):872-79.
- Rathinam VAK, *et al.*, 2019. Innate immunity to intracellular LPS. *Nat. Immunol.* 20(5): 527-533.
- Mather J. & Roberts E. 1998. Contamination: How to avoid it, recognize it, and get rid of it. In: *Introduction to cell and tissue culture: theory and technique.* Chapt 7. p99-9.

# Mycoplasma?

**NEW**

## MycoStrip™ Mycoplasma Detection Kit

*A new way to detect mycoplasma in your cell culture*

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