TABLE OF CONTENTS

Introduction

STING activity

STING agonists

STING signaling inhibitors

STING regulation

Genetic variation in STING

Therapeutic targeting of STING

Conclusion

InvivoGen

EUROPE

5, rue Jean Rodier
F-31400 Toulouse
France
Tel: +33 (0)5 62 71 69 39
Fax: +33 (0)5 62 71 69 30
info.eu@invivogen.com

USA

10515 Vista Sorrento Parkway
San Diego, CA 92121
Tel: +1 888 457 5873
Fax: +1 858 457 5843
info@invivogen.com

ASIA

Unit 106, 1F,
15W Phase 3, Hong Kong Science Park,
Pak Shek Kok,
Hong Kong
Tel: +852 3622 3480
Fax: +852 3622 3483
info.hk@invivogen.com

www.invivogen.com

Follow the path to STING

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

Figure 1: The STING signaling pathway

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
Table 1: Reporter cell lines related to cGAS/STING signaling

CELL LINE	PRODUCTS	DESCRIPTION	UNIT SIZE	CAT.CODE
B16	B16-Blue™ ISG Cells	IRF-SEAP reporter mouse melanoma cells	3-7 x 10 ⁶ cells	bb-ifnabg
	B16-Blue™ ISG-KO-STING Cells	IRF-SEAP reporter STING knockout cells	3-7 x 10 ⁶ cells	bb-kostg
HEK293	HEK-Blue™ ISG Cells	IRF-SEAP reporter human embryonic kidney cells	3-7 x 10 ⁶ cells	hkb-isg
	HEK-Blue™ ISG-KO-STING Cells	IRF-SEAP reporter STING knockout cells	3-7 x 10 ⁶ cells	hkb-kostg
HEK293T	293T-Dual™ hSTING-A162 Cells	IRF-SEAP and IFN- β -Lucia reporter cells with A162 human STING	3-7 x 10 ⁶ cells	293d-a162
	293T-Dual™ hSTING-H232 Cells	IRF-SEAP and IFN- β -Lucia reporter cells with H232 human STING	3-7 x 10 ⁶ cells	293d-h232
	293T-Dual™ hSTING-R232 Cells	IRF-SEAP and IFN- β -Lucia reporter cells with R232 human STING	3-7 x 10 ⁶ cells	293d-r232
	293T-Dual™ mSTING Cells	IRF-SEAP and IFN- β -Lucia reporter cells with murine STING	3-7 x 10 ⁶ cells	293d-mstg
RAW 264.7	RAW-Lucia™ ISG Cells	IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	rawl-isg
	RAW-Lucia™ ISG-KO-cGAS Cells	IRF-Lucia reporter cGAS knockout cells	3-7 x 10 ⁶ cells	rawl-kocgas
	RAW-Lucia™ ISG-KO-IRF3 Cells	IRF-Lucia reporter IRF3 knockout cells	3-7 x 10 ⁶ cells	rawl-koirf3
	RAW-Lucia™ ISG-KO-STING Cells	IRF-Lucia reporter STING knockout cells	3-7 x 10 ⁶ cells	rawl-kostg
THP-1	THP1-Dual™ Cells	NF- κ B-SEAP and IRF-Lucia reporter human monocyte cells	3-7 x 10 ⁶ cells	thpd-nfs
	THP1-Dual™ KO-cGAS Cells	NF- κ B-SEAP and IRF-Lucia reporter cGAS knockout cells	3-7 x 10 ⁶ cells	thpd-kocgas
	THP1-Dual™ KO-STING Cells	NF- κ B-SEAP and IRF-Lucia reporter STING knockout cells	3-7 x 10 ⁶ cells	thpd-kostg
	THP1-Dual™ KI-hSTING-A162 Cells	NF- κ B-SEAP and IRF-Lucia reporter A162 human STING knockin cells	3-7 x 10 ⁶ cells	thpd-a162
	THP1-Dual™ KI-hSTING-H232 Cells	NF- κ B-SEAP and IRF-Lucia reporter H232 human STING knockin cells	3-7 x 10 ⁶ cells	thpd-h232
	THP1-Dual™ KI-hSTING-M155 Cells	NF- κ B-SEAP and IRF-Lucia reporter M155 human STING knockin cells	3-7 x 10 ⁶ cells	thpd-m155
	THP1-Dual™ KI-hSTING-R232 Cells	NF- κ B-SEAP and IRF-Lucia reporter R232 human STING knockin cells	3-7 x 10 ⁶ cells	thpd-r232
	THP1-Dual™ KI-hSTING-S154 Cells	NF- κ B-SEAP and IRF-Lucia reporter S154 human STING knockin cells	3-7 x 10 ⁶ cells	thpd-s154

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
Figure 2: Induction of the interferon regulatory factor pathway by various STING ligands in THP1-Dual™ cells. IRF induction was determined by measuring the relative light units (RLUs) in a luminometer using QUANTI-Luc™, a Lucia luciferase detection reagent. The IRF induction of each ligand is expressed relative to that of hIFN- β at 1×10^4 U/ml (taken as 100%).

Table 2: Cyclic dinucleotides and DMXAA

CATEGORY	PRODUCTS	DESCRIPTION	UNIT SIZE	CAT.CODE
cGAMP	3'3'-cGAMP	Cyclic [G(3',5')pA(3',5')p]	500 μ g	tlrl-nacga
	3'3'-cGAMP VacciGrade™	Preclinical grade of cyclic [G(3',5')pA(3',5')p]	1 mg	vac-nacga
	2'3'-cGAMP	Cyclic [G(2',5')pA(3',5')p]	500 μ g	tlrl-nacga23
	2'3'-cGAMP VacciGrade™	Preclinical grade of cyclic [G(2',5')pA(3',5')p]	1 mg	vac-nacga23
	2'3'-cGAM(PS)₂ (Rp/Sp)	Bisphosphorothioate analog of 2'3'-cGAMP	250 μ g	tlrl-nacga2srs
c-di-GMP	c-di-GMP	Cyclic [G(3',5')pG(3',5')p]	1 mg	tlrl-nacdg
	c-di-GMP VacciGrade™	Preclinical grade of cyclic [G(3',5')pG(3',5')p]	1 mg	vac-nacdg
	2'3'-c-di-GMP	Analog of c-di-GMP	500 μ g	tlrl-nacdg23
c-di-AMP	c-di-AMP	Cyclic [A(3',5')pA(3',5')p]	1 mg	tlrl-nacda
	c-di-AMP VacciGrade™	Preclinical grade of cyclic [A(3',5')pA(3',5')p]	1 mg	vac-nacda
	2'3'-c-di-AMP	Analog of c-di-AMP	500 μ g	tlrl-nacda23
	2'3'-c-di-AM(PS)₂ (Rp,Rp)	Bisphosphorothioate analog of 2'3'-c-di-AMP	100 μ g	tlrl-nacda2r-01
	2'3'-c-di-AM(PS)₂ (Rp,Rp) VacciGrade™	Preclinical grade of bisphosphorothioate analog of 2'3'-c-di-AMP	500 μ g	vac-nacda2r
c-AIMP	cAIMP	Cyclic [A(3',5')pI(3',5')p]	500 μ g	tlrl-nacai
	cAIMP Difluor	Difluor cyclic [A(3',5')pI(3',5')p]	250 μ g	tlrl-nacaidf
	cAIM(PS)₂ Difluor (Rp/Sp)	Difluor and bisphosphorothioate analog of cAIMP	100 μ g	tlrl-nacairs
Non-CDN	DMXAA	5,6-dimethyl-xanthenone-4-acetic acid	5 mg	tlrl-dmx

Inhibition of cGAS/STING signaling

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Table 3: Synthetic inhibitors of the cGAS/STING pathway

PRODUCTS	DESCRIPTION	UNIT SIZE	CAT.CODE
Amlexanox	TBK1/IKKε inhibitor	50 mg	inh-amx
Bay 11-7082	IκB-α inhibitor - NLRP3 inflammasome inhibitor	10 mg	tlrl-b82
BX795	TBK1/IKKε inhibitor - TLR signaling inhibitor - RLR inhibitor	5 mg	tlrl-bx7
Celastrol	NF-κB inhibitor	1 mg	ant-cls
MG-132	26S proteasome inhibitor - Autophagy activator	5 mg	tlrl-mg132
ODN TTAGGG (A151)	cGAS inhibitor - TLR9 inhibitor - AIM2 inhibitor	1 mg	tlrl-ttag151-1
SB202190	MAP kinase inhibitor - Autophagy inducer	5 mg	tlrl-sb90

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
Figure 3: Regulation of STING signaling.

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* All plasmids are provided as 20 μg of lyophilized DNA.


Figure 4: Schematic of human STING genetic variants and their functional effect or disease association.

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
Figure 5: Central role of STING in sensing nucleic acids and in inflammation.

Abbreviations

CDNs: cyclic dinucleotides
cGAS: cyclic GMP-AMP synthase
ER: endoplasmic reticulum
IFNs: interferons
IRF3: interferon regulatory factor 3
ISGs: interferon-stimulated genes
ISRE: interferon-responsive element
NF-κB: nuclear factor kappa light-chain-enhancer of activated B cells
PRR: pattern recognition receptor
SAVI: STING-Associated Vasculopathy with Onset in Infancy
SNP: single nucleotide polymorphism
STING: stimulator of interferon genes
TBK1: TANK-binding-kinase-1
TNF-α: tumor necrosis factor alpha

References

1. Barbalat R. et al., 2011. Nucleic Acid Recognition by the Innate Immune System. Annual Review of Immunology 29: 185-214.
2. Ishikawa H., and Barber G.N. 2008. STING is an endoplasmic reticulum adaptor that facilitates innate immune signalling. Nature 455: 674-8.
3. Zhong B. et al., 2008. The Adaptor Protein MITA Links Virus-Sensing Receptors to IRF3 Transcription Factor Activation. Immunity 29: 538-50.
4. Burdette D.L. et al., 2011. STING is a direct innate immune sensor of cyclic di-GMP. Nature 478: 515-8.
5. Sun L. et al., 2013. Cyclic GMP-AMP Synthase Is a Cytosolic DNA Sensor That Activates the Type I Interferon Pathway. Science 339: 786-91.
6. Wu J. et al., 2013. Cyclic GMP-AMP Is an Endogenous Second Messenger in Innate Immune Signaling by Cytosolic DNA. Science 339: 826-30.
7. Ishikawa H. et al., 2009. STING regulates intracellular DNA-mediated, type I interferon-dependent innate immunity. Nature 461: 788-92.
8. Abe T. and Barber G.N., 2014. Cytosolic-DNA-Mediated, STING-Dependent Proinflammatory Gene Induction Necessitates Canonical NF-κB Activation through TBK1. Journal of Virology 88:5328-41.
9. Karaolis D.K.R. et al., 2005. c-di-GMP (3'-5'-Cyclic Diguanylic Acid) Inhibits Staphylococcus aureus Cell-Cell Interactions and Biofilm Formation. Antimicrobial Agents and Chemotherapy 49: 1029-38.
10. Ogunniyi A.D. et al., 2008. c-di-GMP is an effective immunomodulator and vaccine adjuvant against pneumococcal infection. Vaccine 26: 4676-85.
11. Amikam D. et al., 1995. The novel cyclic dinucleotide 3'-5' cyclic diguanylic acid binds to p21ras and enhances DNA synthesis but not cell replication in the Molt 4 cell line. Biochemical Journal 311: 921-7.
12. Karaolis D.K.R. et al., 2005. 3',5'-Cyclic diguanylic acid (c-di-GMP) inhibits basal and growth factor-stimulated human colon cancer cell proliferation. Biochemical and Biophysical Research Communications 329: 40-5.
13. Steinberger O. et al., 1999. Elevated expression of the CD4 receptor and cell cycle arrest are induced in Jurkat cells by treatment with the novel cyclic dinucleotide 3',5'-cyclic diguanylic acid. FEBS Letters 444: 125-9.
14. Lioux T. et al., 2016. Design, Synthesis, and Biological Evaluation of Novel Cyclic Adenosine-Inosine Monophosphate (cAIMP) Analogs That Activate Stimulator of Interferon Genes (STING). Journal of Medicinal Chemistry 59: 10253-67.
15. Parvatiyar K. et al., 2012. The helicase DDX41 recognizes the bacterial secondary messengers cyclic di-GMP and cyclic di-AMP to activate a type I interferon immune response. Nat Immunol 13: 1155-61.
16. Omura H. et al., 2016. Structural and Functional Analysis of DDX41: a bispecific immune receptor for DNA and cyclic dinucleotide. Scientific Reports 6: 34756.
17. McFarland A.P. et al., 2017. Sensing of Bacterial Cyclic Dinucleotides by the Oxidoreductase RECON Promotes NF-κB Activation and Shapes a Proinflammatory Antibacterial State. Immunity 46: 433-45.
18. Ching L.M. et al., 1991. Haematological effects in mice of the antitumour agents xanthenone-4-acetic acid, 5,6-dimethyl-xanthenone-4-acetic acid [correction of 5,6-methyl-] and flavone acetic acid. Cancer chemotherapy and pharmacology 28: 414-9.
19. Lara Jr P.N. et al., 2011. Randomized Phase III Placebo-Controlled Trial of Carboplatin and Paclitaxel With or Without the Vascular Disrupting Agent Vandimezan (ASA404) in Advanced Non-Small-Cell Lung Cancer. Journal of Clinical Oncology 29: 2965-71.
20. Prantner D. et al., 2012. 5,6-Dimethylxanthenone-4-acetic Acid (DMXAA) Activates Stimulator of Interferon Gene (STING)-dependent Innate Immune Pathways and Is Regulated by Mitochondrial Membrane Potential. Journal of Biological Chemistry 287: 39776-88.
21. Conlon J. et al., 2013. Mouse, but not Human STING, Binds and Signals in Response to the Vascular Disrupting Agent 5,6-Dimethylxanthenone-4-Acetic Acid. The Journal of Immunology 190: 5216-25.
22. Liu Y. et al., 2014. Activated STING in a Vascular and Pulmonary Syndrome. New England Journal of Medicine 371: 507-18.
23. Corrales L. et al., 2015. Direct Activation of STING in the Tumor Microenvironment Leads to Potent and Systemic Tumor Regression and Immunity. Cell Reports 11: 1018-30.
24. Dey B. et al., 2015. A bacterial cyclic dinucleotide activates the cytosolic surveillance pathway and mediates innate resistance to tuberculosis. Nature medicine 21: 401-6.
25. Kim S. et al., 2013. Anticancer Flavonoids Are Mouse-Selective STING Agonists. ACS Chemical Biology 8: 1396-401.
26. Dey R.J. et al., 2017. Inhibition of innate immune cytosolic surveillance by an M. tuberculosis phosphodiesterase. Nat Chem Biol 13: 210-7.
27. Burdette D.L., and Vance R.E. 2013. STING and the innate immune response to nucleic acids in the cytosol. Nat Immunol 14: 19-26.
28. Lau L. et al., 2015. DNA tumor virus oncogenes antagonize the cGAS-STING DNA-sensing pathway. Science 350: 568-71.
29. Aguirre S. et al., 2017. Dengue virus NS2B protein targets cGAS for degradation and prevents mitochondrial DNA sensing during infection. Nature microbiology 2: 17037.
30. Liu Y. et al., 2015. Hepatitis B Virus Polymerase Disrupts K63-Linked Ubiquitination of STING To Block Innate Cytosolic DNA-Sensing Pathways. Journal of Virology 89: 2287-300.
31. Christensen M.H. et al., 2016. HSV-1 ICP27 targets the TBK1-activated STING signalsome to inhibit virus-induced type I IFN expression. The EMBO Journal 35: 1385-99.
32. Fu Y.Z. et al., 2017. Human Cytomegalovirus Tegument Protein UL82 Inhibits STING-Mediated Signaling to Evade Antiviral Immunity. Cell Host & Microbe 21: 231-43.
33. Holm C.K. et al., 2016. Influenza A virus targets a cGAS-independent STING pathway that controls enveloped RNA viruses. Nat Commun 7: 10680.
34. Dobbs N. et al., 2015. STING Activation by Translocation from the ER Is Associated with Infection and Autoinflammatory Disease. Cell Host & Microbe 18: 157-68.
35. Steinhagen F. et al., 2017. Suppressive oligodeoxynucleotides containing TTAGGG motifs inhibit cGAS activation in human monocytes. Eur J Immunol. Dec 7 [Epub ahead of print].
36. Mukai K. et al., 2016. Activation of STING requires palmitoylation at the Golgi. Nature Communications 7: 11932.
37. Pokatayev V. et al., 2016. RNase H2 catalytic core Aicardi-Goutières syndrome-related mutant invokes cGAS-STING innate immune-sensing pathway in mice. The Journal of Experimental Medicine 213: 329-

36. **Chen Y. et al., 2017.** P38 Inhibition Provides Anti-DNA Virus Immunity by Regulation of USP21 Phosphorylation and STING Activation. *J Exp Med.* 214: 991-1010. **39. Dutta D. et al., 2015.** BRCA1 Regulates IFI16 Mediated Nuclear Innate Sensing of Herpes Viral DNA and Subsequent Induction of the Innate Inflammasome and Interferon- β Responses. *PLoS Pathogens* 11: e1005030. **40. Diner B.A. et al., 2015.** The functional interactome of PYHIN immune regulators reveals IFIX is a sensor of viral DNA. *Molecular Systems Biology* 11: 787. **41. Zhang Z. et al., 2011.** The helicase DDX41 senses intracellular DNA mediated by the adaptor STING in dendritic cells. *Nature immunology* 12: 959-65. **42. Liu C. et al., 2016.** AIM2 inhibits autophagy and IFN- β production during *M. bovis* infection. *Oncotarget* 7: 46972-87. **43. Malik P. et al., 2014.** NET23/STING Promotes Chromatin Compaction from the Nuclear Envelope. *PLoS ONE* 9: e111851. **44. Parkes E.E. et al., 2017.** Activation of STING-Dependent Innate Immune Signaling By S-Phase-Specific DNA Damage in Breast Cancer. *JNCI: Journal of the National Cancer Institute* 109: 199. **45. Kondo T. et al., 2013.** DNA damage sensor MRE11 recognizes cytosolic double-stranded DNA and induces type I interferon by regulating STING trafficking. *Proceedings of the National Academy of Sciences* 110: 2969-74. **46. Bhattacharya S. et al., 2017.** RAD51 interconnects between DNA replication, DNA repair and immunity. *Nucleic Acids Research* 45: 4590-605. **47. Pépin G. et al., 2016.** Cre-dependent DNA recombination activates a STING-dependent innate immune response. *Nucleic Acids Research* 44: 5356-64. **48. Wang Y. et al., 2017.** Inflammasome Activation Triggers Caspase-1-Mediated Cleavage of cGAS to Regulate Responses to DNA Virus Infection. *Immunity* 46: 393-404. **49. Chen Q. et al., 2016.** Regulation and function of the cGAS-STING pathway of cytosolic DNA sensing. *Nat Immunol* 17: 1142-9. **50. Maringer K., and Fernandez-Sesma A. 2014.** Message in a bottle: lessons learned from antagonism of STING signalling during RNA virus infection. *Cytokine & Growth Factor Reviews* 25: 669-79. **51. Zevini A. et al., 2017.** Crosstalk between Cytosolic RIG-I and STING Sensing Pathways. *Trends in Immunology* 38: 194-205. **52. Suspène R. et al., 2017.** Self-cytoplasmic DNA upregulates the mutator enzyme APOBEC3A leading to chromosomal DNA damage. *Nucleic Acids Research* 45: 3231-41. **53. Liu Y. et al., 2017.** RIGulation of STING expression: at the crossroads of viral RNA and DNA sensing pathways. *Inflammation and cell signaling* 4: e1491. **54. Goulet M.-L. et al., 2013.** Systems Analysis of a RIG-I Agonist Inducing Broad Spectrum Inhibition of Virus Infectivity. *PLOS Pathogens* 9: e1003298. **55. Härtlova A. et al., 2015.** DNA Damage Primes the Type I Interferon System via the Cytosolic DNA Sensor STING to Promote Anti-Microbial Innate Immunity. *Immunity* 42: 332-43. **56. Yu X. et al., 2016.** Cross-Regulation of Two Type I Interferon Signaling Pathways in Plasmacytoid Dendritic Cells Controls Anti-malaria Immunity and Host Mortality. *Immunity* 45: 1093-107. **57. Wang X. et al., 2016.** STING Requires the Adaptor TRIF to Trigger Innate Immune Responses to Microbial Infection. *Cell Host & Microbe* 20: 329-41. **58. Chen M. et al., 2016.** TRIM14 Inhibits cGAS Degradation Mediated by Selective Autophagy Receptor p62 to Promote Innate Immune Responses. *Molecular Cell* 64: 105-19. **59. Moretti J. et al., 2017.** STING senses microbial viability to orchestrate stress-mediated autophagy of the endoplasmic reticulum. *Cell* 171:1-15. **60. Shibutani S.T. et al., 2015.** Autophagy and autophagy-related proteins in the immune system. *Nat Immunol* 16: 1014-24. **61. Watson R.O. et al., 2015.** The cytosolic sensor cGAS detects Mycobacterium tuberculosis DNA to induce type I interferons and activate autophagy. *Cell host & microbe* 17: 811-9. **62. Liang Q. et al., 2014.** Crosstalk between the cGAS DNA Sensor and Beclin-1 Autophagy Protein Shapes Innate Antimicrobial Immune Responses. *Cell Host & Microbe* 15: 228-38. **63. Saitoh T. et al., 2009.** Atg9a controls dsDNA-driven dynamic translocation of STING and the innate immune response. *Proceedings of the National Academy of Sciences* 106: 20842-6. **64. Petrasek J. et al., 2013.** STING-IRF3 pathway links endoplasmic reticulum stress with hepatocyte apoptosis in early alcoholic liver disease. *Proceedings of the National Academy of Sciences* 110: 16544-9. **65. Iracheta-Vellve A. et al., 2016.** Endoplasmic Reticulum Stress-induced Hepatocellular Death Pathways Mediate Liver Injury and Fibrosis via Stimulator of Interferon Genes. *Journal of Biological Chemistry* 291: 26794-805. **66. White Michael J. et al., 2014.** Apoptotic Caspases Suppress mtDNA-Induced STING-Mediated Type I IFN Production. *Cell* 159: 1549-62. **67. Tang C.-H.A. et al., 2016.** Agonist-Mediated Activation of STING Induces Apoptosis in Malignant B Cells. *Cancer Research* 76: 2137-52. **68. Cui Y. et al., 2016.** Mycobacterium bovis Induces Endoplasmic Reticulum Stress Mediated-Apoptosis by Activating IRF3 in a Murine Macrophage Cell Line. *Frontiers in Cellular and Infection Microbiology* 6: 182. **69. Hu M.-M. et al., 2016.** Sumoylation Promotes the Stability of the DNA Sensor cGAS and the Adaptor STING to Regulate the Kinetics of Response to DNA Virus. *Immunity* 45: 555-69. **70. Chiang C., and Gack M.U. 2017.** Post-translational Control of Intracellular Pathogen Sensing Pathways. *Trends in Immunology* 38: 39-52. **71. Wang P-H. et al., 2018.** A novel transcript isoform of STING that sequesters cGAMP and dominantly inhibits innate nucleic acid sensing. *Nucleic Acids Res. Mar 14* [Epub ahead of print]. **72. Yi G. et al., 2013.** Single Nucleotide Polymorphisms of Human STING Can Affect Innate Immune Response to Cyclic Dinucleotides. *PLOS ONE* 8: e77846. **73. Patel S. et al., 2017.** The Common R71H-G230A-R293Q Human TMEM173 Is a Null Allele. *The Journal of Immunology* 198: 776-87. **74. Chen H. et al., 2014.** An Alternative Splicing Isoform of MITA Antagonizes MITA-Mediated Induction of Type I IFNs. *The Journal of Immunology* 192: 1162-70. **75. Jeremiah N. et al., 2014.** Inherited STING-activating mutation underlies a familial inflammatory syndrome with lupus-like manifestations. *The Journal of Clinical Investigation* 124: 5516-20. **76. Munoz J. et al., 2015.** Stimulator of interferon genes-associated vasculopathy with onset in infancy : A mimic of childhood granulomatosis with polyangiitis. *JAMA Dermatology* 151: 872-7. **77. Warner J.D. et al., 2017.** STING-associated vasculopathy develops independently of IRF3 in mice. *J Exp Med.* 214: 3279-92. **78. Melki I. et al., 2017.** Disease-associated mutations identify a novel region in human STING necessary for the control of type I interferon signaling. *Journal of Allergy and Clinical Immunology* 140: 543-52.e5. **79. König N. et al., 2017.** Familial chilblain lupus due to a gain-of-function mutation in STING. *Annals of the Rheumatic Diseases* 76: 468-72. **80. Dubensky T.W. et al., 2013.** Rationale, progress and development of vaccines utilizing STING-activating cyclic dinucleotide adjuvants. *Therapeutic Advances in Vaccines* 1: 131-43. **81. Ma Z., and Damania B. 2016.** The cGAS-STING Defense Pathway and Its Counteraction by Viruses. *Cell Host & Microbe* 19: 150-8. **82. Gentili M. et al., 2015.** Transmission of innate immune signaling by packaging of cGAMP in viral particles. *Science* 349: 1232-6. **83. Ablasser A. et al., 2013.** cGAS produces a 2'-5'-linked cyclic dinucleotide second messenger that activates STING. *Nature* 498: 380-4. **84. Guo H. et al., 2016.** NLRX1 sequesters STING to negatively regulate the interferon response, thereby facilitating the replication of HIV-1 and DNA viruses. *Cell Host and Microbes.* 19: 515-28. **85. Mitzel D.N. et al., 2014.** Age-Enhanced Endoplasmic Reticulum Stress Contributes to Increased Atg9A Inhibition of STING-Mediated IFN- β Production during *Streptococcus pneumoniae* Infection. *The Journal of Immunology* 192: 4273-83. **86. Hansen C. et al., 2014.** Listeria monocytogenes induces IFN β expression through an IFI16-, cGAS- and STING-dependent pathway. *The EMBO Journal* 33: 1654-66. **87. Archer K.A. et al., 2014.** STING-Dependent Type I IFN Production Inhibits Cell-Mediated Immunity to Listeria monocytogenes. *PLOS Pathogens* 10: e1003861. **88. Chandra D. et al., 2014.** STING Ligand c-di-GMP Improves Cancer Vaccination against Metastatic Breast Cancer. *Cancer Immunology Research* 2: 901-10. **89. Xia T. et al., 2016.** Deregulation of STING Signaling in Colorectal Carcinoma Constrains DNA Damage Responses and Correlates With Tumorigenesis. *Cell Reports* 14: 282-97. **90. Demaria O. et al., 2015.** STING activation of tumor endothelial cells initiates spontaneous and therapeutic antitumor immunity. *Proceedings of the National Academy of Sciences* 112: 15408-13. **91. Deng L. et al., 2014.** STING-Dependent Cytosolic DNA Sensing Promotes Radiation-Induced Type I Interferon-Dependent Antitumor Immunity in Immunogenic Tumors. *Immunity* 41: 843-52. **92. <https://clinicaltrials.gov/ct2/show/NCT02675439?term=adu-s100&rank=1>** **93. Ahn J., Konno H., and Barber G.N. 2015.** Diverse roles of STING-dependent signaling on the development of cancer. *Oncogene* 34: 5302-8. **94. Zhu Q. et al., 2014.** Cutting Edge: STING Mediates Protection against Colorectal Tumorigenesis by Governing the Magnitude of Intestinal Inflammation. *The Journal of Immunology* 193: 4779-82. **95. Lemos H. et al., 2016.** STING Promotes the Growth of Tumors Characterized by Low Antigenicity via IDO Activation. *Cancer Research* 76: 2076-81. **96. Starokadomskyy P. et al., 2016.** DNA polymerase- α regulates the activation of type I interferons through cytosolic RNA:DNA synthesis. *Nat Immunol* 17: 495-504. **97. An J. et al., 2017.** Expression of Cyclic GMP-AMP Synthase in Patients With Systemic Lupus Erythematosus. *Arthritis & Rheumatology* 69: 800-7. **98. York Autumn G. et al., 2015.** Limiting Cholesterol Biosynthetic Flux Spontaneously Engages Type I IFN Signaling. *Cell* 163: 1716-29.

InvivoGen
intofocus

For more information
please visit : www.invivogen.com

STING: STIMULATOR OF INTERFERON GENES

EUROPE

5, rue Jean Rodier
F-31400 Toulouse
France
Tel: +33 (0)5 62 71 69 39
Fax: +33 (0)5 62 71 69 30
info.eu@invivogen.com

USA

10515 Vista Sorrento Parkway
San Diego, CA 92121
Tel: +1 888 457 5873
Fax: +1 858 457 5843
info@invivogen.com

ASIA

Unit 106, 1F, 15W Phase 3,
Hong Kong Science Park,
Pak Shek Kok,
Hong Kong
Tel: +852 3622 3480
Fax: +852 3622 3483
info.hk@invivogen.com