INVIVOGEN'S PRODUCTS FOCUS

STING LIGANDS

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cAIM(PS)2 Difluor (Rp/Sp) (CL656) **RECOMMENDED!**

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XANTHENONE ANALOG

DMXAA (also known as Vadimezan or ASA404)

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STING variants differentially respond to CDN


Fig 1: The IRF pathway (A) and NF- κ B pathway (B) induction in THP1-Dual™ KI-hSTING-H232 cells in response to the STING ligands CDNs has been assessed after 24h incubation. In contrast to THP1-Dual™ KI-hSTING-R232 cells expressing the “wild-type” R232 human STING variant, THP1-Dual™ KI-hSTING-H232 cells do not respond to bacterial CDNs such as 3'3'-cGAMP and display only a weak IRF induction in response to mammalian CDNs such as 2'3'-cGAMP, while 3'3'-cGAMP Fluorinated or 2'3'-cGAMP(PS)₂ (Rp/Sp) induce a moderate IRF induction (red arrows). In THP1-Dual™ KI-hSTING-H232 cells, no NF- κ B induction was observed following a 24-hour incubation with various CDNs.

IFN induction improved by fluorinated modification


Fig 2: Induction of the interferon regulatory factor (IRF) pathway by various STING ligands in THP1-Dual™ cells. Cells were stimulated for 24 hours with the STING ligands as shown (all at 10 μ g/ml) and assessed IRF induction by measuring the relative light units (RLUs) in a luminometer using QUANTI-Luc™.

IFN induction improved by modification with cAIMPs


Fig 3: Induction of the interferon regulatory factor (IRF) pathway by various STING ligands in THP1-Dual™ cells. Cells were stimulated for 24 hours with CDN at 10 μ g/ml. IRF induction was determined by measuring the relative light units (RLUs) in a luminometer using QUANTI-Luc™.

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Check out InvivoGen's Infocus Newsletter to learn more about STING

InvivoGen
infocus

TABLE OF CONTENTS

- Introduction
- STING activity
- STING-specific inhibitors
- STING regulation
- Genetic variation in STING
- Therapeutic targeting of STING
- Conclusion

InvivoGen
2015-2016
The InvivoGen Infocus Newsletter
is available at www.invivogen.com/infocus

Use STING, the InvivoGen Center
for Research and Development
at www.invivogen.com/sting

Follow the path to STING

STING (Stimulator of Interferon Genes) has become a focal point in immunology research as well as a target in drug discovery. As a signaling hub in innate immunity, STING is a pattern recognition receptor (PRR) of pericellular lipopeptides in endoplasmic reticulum (ER) and cytoplasm. InvivoGen offers a growing variety of products to help you explore STING, its signaling partners, cytokine induction activity and therapeutic potential.

Introduction

Entry of self or foreign nucleic acids into the cytoplasm can signal various problems, including pathogen infection or tissue damage, mitochondrial or nuclear damage, and the presence of tumors. Until the discovery of STING in 2006, detection of nucleic acids in the cytoplasm was thought to be the domain of the innate immune system. However, STING has been largely ignored by the immune system. In 2013, a study of STING activity in the nucleus and its role in the regulation of type I interferon (IFN) responses and the control of cancer via iNOS/IFN- γ was presented to be an essential molecule which integrates signaling between signaling of both the ER and cytoplasmic STING (cSTING) and its downstream signaling pathway. STING is a novel sensor of cytosolic nucleic acids, which is involved in the regulation of innate immunity. In this review, we discuss the role of STING in the regulation of innate immunity and its potential as a therapeutic target.

InvivoGen STING

InvivoGen offers a growing variety of products to help you explore STING, its signaling partners, cytokine induction activity and therapeutic potential.

STING activity

STING is a pattern recognition receptor (PRR) of pericellular lipopeptides in endoplasmic reticulum (ER) and cytoplasm. InvivoGen offers a growing variety of products to help you explore STING, its signaling partners, cytokine induction activity and therapeutic potential.

STING-specific inhibitors

STING is a pattern recognition receptor (PRR) of pericellular lipopeptides in endoplasmic reticulum (ER) and cytoplasm. InvivoGen offers a growing variety of products to help you explore STING, its signaling partners, cytokine induction activity and therapeutic potential.

STING regulation

STING is a pattern recognition receptor (PRR) of pericellular lipopeptides in endoplasmic reticulum (ER) and cytoplasm. InvivoGen offers a growing variety of products to help you explore STING, its signaling partners, cytokine induction activity and therapeutic potential.

Genetic variation in STING

STING is a pattern recognition receptor (PRR) of pericellular lipopeptides in endoplasmic reticulum (ER) and cytoplasm. InvivoGen offers a growing variety of products to help you explore STING, its signaling partners, cytokine induction activity and therapeutic potential.

Therapeutic targeting of STING

STING is a pattern recognition receptor (PRR) of pericellular lipopeptides in endoplasmic reticulum (ER) and cytoplasm. InvivoGen offers a growing variety of products to help you explore STING, its signaling partners, cytokine induction activity and therapeutic potential.

Conclusion

STING is a pattern recognition receptor (PRR) of pericellular lipopeptides in endoplasmic reticulum (ER) and cytoplasm. InvivoGen offers a growing variety of products to help you explore STING, its signaling partners, cytokine induction activity and therapeutic potential.

Figure 1: The STING signaling pathway

Visit

<https://www.invivogen.com/sting-products>
for product details

PRODUCTS

STING Ligands

PRODUCT	DESCRIPTION	QTY	CAT. CODE
Natural CDNs			
2'3'-cGAMP	Cyclic [G(2',5')pA(3',5')p]	200 µg 500 µg 1 mg 5 mg	tlrl-nacga23-02 tlrl-nacga23 tlrl-nacga23-1 tlrl-nacga23-5
2'3'-cGAMP VacciGrade™	Preclinical grade of 2'3'-cGAMP	1 mg	vac-nacga23
3'3'-cGAMP	Cyclic [G(3',5')pA(3',5')p]	500 µg 1 mg 5x 0.5 mg	tlrl-nacga tlrl-nacga-1 tlrl-nacga-2.5
c-di-AMP	Cyclic diadenylate monophosphate	1 mg 5 x 1 mg	tlrl-nacda tlrl-nacda-5
c-di-AMP VacciGrade™	Preclinical grade of c-di-AMP	1 mg	vac-nacda
c-di-GMP	Cyclic diguanylate monophosphate	1 mg 5 x 1 mg	tlrl-nacdg tlrl-nacdg-5
c-di-GMP VacciGrade™	Preclinical grade of c-di-GMP	1 mg	vac-nacdg
cGAMP-derived CDNs			
2'3'-cGAM(PS)2 (Rp/Sp)	Bisphosphorothioate analog of 2'3'-cGAMP, Rp,Sp-isomers	250 µg	tlrl-nacga2srs
3'3'-cGAMP Fluorinated	Difluoro 3'3'-cGAMP	2 x 100 µg	tlrl-nacgaf-2
c-di-AMP-derived CDNs			
c-di-AMP Fluorinated	Difluoro 3'3'-c-di-AMP	2 x 100 µg	tlrl-nacdaf-2
2'3'-c-di-AM(PS)2 (Rp,Rp) (Identical to ADU-S100)	Bisphosphorothioate analog of 2'3'-c-di-AMP, Rp isomers	100 µg 500 µg	tlrl-nacda2r-01 tlrl-nacda2r
2'3'-c-di-AM(PS)2 (Rp,Rp) VacciGrade™	Preclinical grade of 2'3'-c-di-AM(PS)2 (Rp,Rp)	500 µg	vac-nacda2r
c-di-GMP-derived CDNs			
c-di-GMP Fluorinated	Difluoro 3'3'-c-di-GMP	2 x 100 µg	tlrl-nacdgf-2
cAIM-derived CDNs			
cAIMP (CL592)	Cyclic (adenine monophosphate- inosine monophosphate)	500 µg	tlrl-nacai
cAIMP Difluor (CL614)	cAIMP difluorinated	250 µg	tlrl-nacaidf
cAIM(PS)2 Difluor (Rp/Sp) (CL656)	cAIMP bisphosphorothioate and difluorinated, Rp,Sp-isomers	2 x 100 µg	tlrl-nacairs-2
Non-CDN			
DMXAA	Murine STING ligand - Xanthenone Analog	5 mg	tlrl-dmx

Reporter cell lines related to cGAS/STING signaling

PRODUCT	DESCRIPTION	QTY	CAT. CODE
THP-1 monocytes			
THP1-Dual™ Cells	NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-nfis
THP1-Dual™ KO-STING Cells	STING knockout NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-kostg
THP1-Dual™ KO-cGAS Cells	cGAS knockout NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-kocgas
THP1-Dual™ KO-IFI16 cells	IFI16 knockout NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-koif16
THP1-Dual™ KO-IRF3 Cells	IRF3 knockout NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-koirf3
THP1-Dual™ KO-TBK1 Cells	TBK1 knockout NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-kotbk
THP1-Dual™ KO-TREX1 Cells	TREX1 knockout NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-kotrex
THP1-Dual™ KI-hSTING-S154 Cells	NF-κB-SEAP and IRF-Lucia reporter S154 human STING knockin cells	3-7 x 10 ⁶ cells	thpd-s154
THP1-Dual™ KI-hSTING-M155 Cells	NF-κB-SEAP and IRF-Lucia reporter M155 human STING knockin cells	3-7 x 10 ⁶ cells	thpd-m155
THP1-Dual™ KI-hSTING-A162 Cells	NF-κB-SEAP and IRF-Lucia reporter A162 human STING knockin cells	3-7 x 10 ⁶ cells	thpd-a162
THP1-Dual™ KI-hSTING-H232 Cells	STING (H232 isoform) knockin NF-κB-SEAP and IRF-Lucia Reporter Cells	3-7 x 10 ⁶ cells	thpd-h232
THP1-Dual™ KI-hSTING-R232 Cells	STING (R232 isoform) knockin NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-r232
THP1-Dual™ KI-mSTING Cells	Murine STING knockin NF-κB-SEAP and IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	thpd-mstg
293T			
293-Dual™ hSTING-A162 Cells	Dual IRF and IFN-β reporter 293 cells expressing A162 isoform of human STING (S162A)	3-7 x 10 ⁶ cells	293d-a162
293-Dual™ hSTING-H232 Cells	Dual IRF and IFN-β reporter 293 cells expressing H232 isoform of human STING (R232H)	3-7 x 10 ⁶ cells	293d-h232
293-Dual™ hSTING-R232 Cells	Dual IRF and IFN-β reporter 293 cells expressing R232 isoform of human STING	3-7 x 10 ⁶ cells	293d-r232
293-Dual™ mSTING Cells	Dual IRF and IFN-β reporter 293 cells expressing murine STING	3-7 x 10 ⁶ cells	293d-mstg
293-Dual™ Null Cells	Dual IRF and IFN-β reporter 293 cells	3-7 x 10 ⁶ cells	293d-null
HEK 293			
HEK-Blue™ ISG Cells	IRF-inducible SEAP reporter cells	3-7 x 10 ⁶ cells	hkb-isg
HEK-Blue™ ISG-KO-STING Cells	STING knockout IRF-SEAP reporter cells	3-7 x 10 ⁶ cells	hkb-kostg
HEK-Blue™ STAT6-hSTING-R232 Cells	Human STING (R232 variant)-dependent STAT6 HEK293 reporter cells	3-7 x 10 ⁶ cells	hkb-st6r232
Murine B16 melanocytes			
B16-Blue™ ISG Cells	IRF-inducible SEAP reporter cells	3-7 x 10 ⁶ cells	bb-ifnabg
B16-Blue™ ISG-KO-STING Cells	STING knockout IRF-SEAP reporter cells	3-7 x 10 ⁶ cells	bb-kostg
Murine RAW 264.7 macrophages			
RAW-Lucia™ ISG Cells	IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	rawl-isg
RAW-Lucia™ ISG-KO-STING Cells	STING knockout IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	rawl-kostg
RAW-Lucia™ ISG-KO-cGAS Cells	cGAS knockout IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	rawl-kocgas
RAW-Lucia™ ISG-KO-IFI16 Cells	IFI16 knockout IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	rawl-koif16
RAW-Lucia™ ISG-KO-IRF3 Cells	IRF3 knockout IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	rawl-koirf3
RAW-Lucia™ ISG-KO-TBK1 Cells	TBK1 knockout IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	rawl-kotbk
RAW-Lucia™ ISG-KO-TREX1 Cells	TREX1 knockout IRF-Lucia reporter cells	3-7 x 10 ⁶ cells	rawl-kotrex

Other ligands related to cGAS/STING signaling

PRODUCT	DESCRIPTION	QTY	CAT. CODE
dsDNA-EC	CDS & TLR9 Agonist - Double-stranded genomic DNA from E. coli K12	200 µg	tlrl-ecdna
G3-YSD	Y-form DNA - cGAS Agonist	200 µg	tlrl-ydna
HSV-60	CDS Agonist	200 µg	tlrl-hsv60n
HSV-60/LyoVec™	Complexed with LyoVec™	100 µg	tlrl-hsv60c
ISD	Interferon stimulatory DNA - CDS Agonist	200 µg	tlrl-isdn
ISD/LyoVec™	Complexed with LyoVec™	100 µg	tlrl-isdc
ODN TTAGGG (A151)	Suppressive oligonucleotide - human preferred cGAS, TLR9 and AIM2 Antagonist	200 µg 1 mg	tlrl-ttag151 tlrl-ttag151-1
Poly(dA:dT)	dsDNA naked - CDS, RIG-I Agonist and AIM2 Inducer	200 µg 1 mg	tlrl-patn tlrl-patn-1
Poly(dA:dT)/LyoVec™	Complexed with LyoVec™	100 µg	tlrl-patc
Poly(dA:dT) Rhodamine	Rhodamine labeled	10 µg	tlrl-patrh
Poly(dG:dC)	dsDNA naked - CDS Agonist	200 µg	tlrl-pgcn
Poly(dG:dC)/LyoVec™	Complexed with LyoVec™	100 µg	tlrl-pgcc
VACV-70	CDS Agonist	200 µg	tlrl-vav70n
VACV-70/LyoVec™	Complexed with LyoVec™	100 µg	tlrl-vav70c

Synthetic inhibitors of cGAS/STING pathway

PRODUCT	DESCRIPTION	QTY	CAT. CODE
BX795	TBK1/IKKε inhibitor	5 mg	tlrl-bx7
G140	Human cGAS inhibitor	2 mg	inh-g140
H-151	Synthetic Indole Derivative - STING Inhibitor	10 mg	inh-h151
RU.521	Murine cGAS inhibitor	2 mg	inh-ru521

Selective antibiotics


PRODUCT	DESCRIPTION	QTY	CAT. CODE
Blasticidin	Selective antibiotic for the bsr, bls, or BSD genes	100 mg (10 x 1 ml)	ant-bl-1
G418 (Geneticin)	Selective antibiotic for the neo gene	1 g (10 x 1 ml)	ant-gn-1
Hygromycin B Gold™	Selective antibiotic for the hph gene	1 g (10 x 1 ml)	ant-hg-1
Puromycin	Selective antibiotic for the pac gene	100 mg (10 x 1 ml)	ant-pr-1
Zeocin™	Selective antibiotic for the Sh ble gene	1 g (10 x 1 ml)	ant-zn-1

REFERENCES

1. Kranzusch, P.J. et al., 2015. Ancient origin of cGAS-STING reveals mechanism of universal 2', 3' cGAMP signaling. *Molecular cell*, 59(6), 891-903.
2. Li L. et al., 2014. Hydrolysis of 2'3'-cGAMP by ENPP1 and design of nonhydrolyzable analogs. *Nat Chem Biol*. 10(12):1043-8.
3. Conlon J. et al., 2013. Mouse, but not human STING, binds and signals in response to the vascular disrupting agent 5,6-dimethylxanthenone-4-acetic acid. *J Immunol*. 190(10):5216-25.
4. Kim S. et al., 2013. Anticancer Flavonoids Are Mouse-Selective STING Agonists. *ACS Chem Biol*. 8(7): 1396-1401.
5. Che X. et al., 2017. Single Mutations Reshape the Structural Correlation Network of the DMXAA-Human STING Complex. *J Phys Chem B*. 2017 Mar 9;121(9):2073-2082.
6. Olagnier, D.P. et al., 2020. Identification of SARS-CoV2-mediated suppression of NRF2 signaling reveals a potent antiviral and anti-inflammatory activity of 4-octyl-itaconate and dimethyl fumarate. [Pre-print]
7. Han, L. et al. 2020. SARS-CoV-2 ORF9b antagonizes type I and III interferons by targeting multiple components of RIG-I/MDA-5-MAVS, TLR3-TRIF, and cGAS-STING signaling pathways. *bioRxiv*.
8. Dey, R. et al., 2017. Inhibition of innate immune cytosolic surveillance by an M. tuberculosis phosphodiesterase. *Nat Chem Biol* 13, 210-217.
9. Mitzel, D.N. et al., 2014. Age-enhanced endoplasmic reticulum stress contributes to increased Atg9A inhibition of STING-mediated IFN- β production during *Streptococcus pneumoniae* infection. *The Journal of Immunology*, 192(9), 4273-4283.
10. Koubou, M.K. et al., 2018. STING agonists enable antiviral cross-talk between human cells and confer protection against genital herpes in mice. *PLoS pathogens*, 14(4), e1006976.
11. Kato, Y. et al., 2018. Apoptosis-derived membrane vesicles drive the cGAS-STING pathway and enhance type I IFN production in systemic lupus erythematosus. *Annals of the rheumatic diseases*, 77(10), 1507-1515.
12. Liu, Y. et al., 2014. Activated STING in a vascular and pulmonary syndrome. *New England Journal of Medicine*, 371(6), 507-518.
13. Tansakul, M. et al., 2020. Deficiency of STING Promotes Collagen-Specific Antibody Production and B Cell Survival in Collagen-Induced Arthritis. *Frontiers in Immunology*, 11, 1101.
14. Pokatayev, V. et al., 2020. Homeostatic regulation of STING protein at the resting state by stabilizer TOLLIP. *Nature immunology*, 21(2), 158-167.
15. Ghaffari, A. et al., 2018. STING agonist therapy in combination with PD-1 immune checkpoint blockade enhances response to carboplatin chemotherapy in high-grade serous ovarian cancer. *British journal of cancer*, 119(4), 440-449.
16. Thomsen, M.K. et al., 2020. The cGAS-STING pathway is a therapeutic target in a preclinical model of hepatocellular carcinoma. *Oncogene*, 39(8), 1652-1664.
17. Wu, J. et al., 2020. Interferon-independent activities of mammalian STING mediate antiviral response and tumor immune evasion. *Immunity*, 53(1), 115-126.
18. Zhou, Y. et al., 2020. Blockade of the phagocytic receptor MerTK on tumor-associated macrophages enhances P2X7R-dependent STING activation by tumor-derived cGAMP. *Immunity*, 52(2), 357-373.
19. Tang, C.H.A. et al., 2016. Agonist-mediated activation of STING induces apoptosis in malignant B cells. *Cancer research*, 76(8), 2137-2152.
20. Xu, N. et al., 2021. STING agonist promotes CAR T cell trafficking and persistence in breast cancer. *Journal of Experimental Medicine*, 218(2).
21. Shi, Y. et al., 2020. Intratumoral accumulation of gut microbiota facilitates CD47-based immunotherapy via STING signaling. *The Journal of experimental medicine*, 217(5), e20192282.
22. Chen, J. et al., 2020. Parenteral immunization with a cyclic guanosine monophosphate-adenosine monophosphate (cGAMP) adjuvanted *Helicobacter pylori* vaccine induces protective immunity against *H. pylori* infection in mice. *Human Vaccines & Immunotherapeutics*, 1-6.
23. Li, X.D. et al., 2013. Pivotal roles of cGAS-cGAMP signaling in antiviral defense and immune adjuvant effects. *Science*, 341(6152), 1390-1394.
24. Luo, J. et al., 2019. Enhancing immune response and heterosubtypic protection ability of inactivated H7N9 vaccine by using STING agonist as a mucosal adjuvant. *Frontiers in immunology*, 10, 2274.
25. Goodwin, T.J. & Huang, L., 2017. Investigation of phosphorylated adjuvants co-encapsulated with a model cancer peptide antigen for the treatment of colorectal cancer and liver metastasis. *Vaccine*, 35(19), 2550-2557.
26. Martin, T.L. et al., 2017. Sublingual targeting of STING with 3' 3'-cGAMP promotes systemic and mucosal immunity against anthrax toxins. *Vaccine*, 35(18), 2511-2519.
27. Muñoz González, F. et al., 2019. The BtaF adhesin is necessary for full virulence during respiratory infection by *Brucella suis* and is a novel immunogen for nasal vaccination against *Brucella* infection. *Frontiers in Immunology*, 10, 1775.
28. Wang, Z. & Celis, E., 2015. STING activator c-di-GMP enhances the anti-tumor effects of peptide vaccines in melanoma-bearing mice. *Cancer Immunology, Immunotherapy*, 64(8), 1057-1066.
29. Ebensen, T. et al., 2011. Bis-(3' 5')-cyclic dimeric adenosine monophosphate: strong Th1/Th2/Th17 promoting mucosal adjuvant. *Vaccine*, 29(32), 5210-5220.
30. Fu, J., et al., 2015. STING agonist formulated cancer vaccines can cure established tumors resistant to PD-1 blockade. *Sci Transl Med*. 7(283), 283ra52.

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