CONTENTS:

EDITORIAL

Importance of mAb murinization for *in vivo* mouse models

Cell-based assays for screening SARS-CoV-2 therapeutics

PRODUCTS

Mouse Anti-Mouse mAbs

- Immune Checkpoints
 - CTLA-4
 - PD-1
 - PD-L1
- Lymphocyte Markers for Depletion
 - CD3
 - CD4
 - CD8
- Tumor-Associated Antigens
 - gp75

COVID-19 Related Cell Lines

- HEK-Blue™ hACE2 Cells
- A549-hACE2 Cells
- A549-hACE2-TMPRSS2 Cells

InvivoGen

ASIA

Unit 106, 1F,
15W Phase 3, Hong Kong Science Park,
Pak Shek Kok,
Hong Kong

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

REFERENCES

1. Almagro, J.C. et al. 2017. Progress and Challenges in the Design and Clinical Development of Antibodies for Cancer Therapy. *Front Immunol* 8, 1751.
2. Kimball, J.A. et al. 1993. OKT3 antibody response study (OARS): a multicenter comparative study. *Transplant Proc* 25, 558-560.
3. Shawler, D.L. et al. 1985. Human immune response to multiple injections of murine monoclonal IgG. *J Immunol* 135, 1530-1535.
4. Adair, J.R. et al. 1994. Humanization of the murine anti-human CD3 monoclonal antibody OKT3. *Hum Antibodies Hybridomas* 5(1-2):41-7.
5. Mall, C. et al. 2016. Repeated PD-1/PD-L1 monoclonal antibody administration induces fatal xenogeneic hypersensitivity reactions in a murine model of breast cancer. *Oncimmunology* 5, e1075114.
6. Murphy, J.T. et al. 2014. Anaphylaxis caused by repetitive doses of a GITR agonist monoclonal antibody in mice. *Blood* 123, 2172-2180.
7. Belmar, N.A. et al. 2017. Murinization and H Chain Isotype Matching of the Anti-GITR Antibody DTA-1 Reduces Immunogenicity-Mediated Anaphylaxis in C57BL/6 Mice. *J Immunol* 198, 4502-4512.

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Murinization for reduced toxicity *in vivo*

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LIMITED RISK OF A HYPERSENSITIVITY REACTION WITH A MURINIZED ANTIBODY

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Recombinant production for increased quality

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


Fig 1: Body temperature comparison between human and murinized Anti-PD-L1 mAb after IP4. BALB/c mice aged 10 weeks were challenged subcutaneously with 2×10^5 CT26 cells. After 8 days, Anti-PD-L1-mIgG1e3 (murine) or Anti-PD-L1-hIgG1NQ (human) was administered intraperitoneally (IP) into the mice in sterile PBS (200 µg/mouse). IP injections were then performed twice a week for 3 weeks. After IP4, rectal temperature measurements were taken every 10 minutes for all groups. Data are presented as mean $\pm$ SEM.

REFERENCES CONTINUED

8. Bradbury, A. et al. 2018. When monoclonal antibodies are not monospecific: Hybridomas frequently express additional functional variable regions. *mAbs*, 10(4), 539–546. 9. Selby, M.J. et al. 2013. Anti-CTLA-4 antibodies of IgG2a isotype enhance antitumor activity through reduction of intratumoral regulatory T cells. *Cancer Immunol Res.* 1(1):32–42.

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
Fig 2: Tumor growth after treatment with a hybridoma-produced or recombinantly-produced Anti-PD-1 mAb. BALB/c mice aged 10 weeks were challenged subcutaneously with 2×10^5 CT26 cells. After 8 days, InvivoGen's Anti-PD-1-mlgG1e3 or a hybridoma-produced Anti-PD-1 mAb was administered intraperitoneally (IP) into the mice in sterile PBS (200 µg/mouse). Following this, IP injections were performed twice a week for 3 weeks (red arrows). Tumor growth was monitored for 60 days.


Fig 3: Tumor growth after treatment with a murine Anti-CTLA4-mlgG2b and Anti-CTLA4-mlgG2a mAb. BALB/c mice aged 10 weeks were challenged subcutaneously with 2×10^5 CT26 cells. After 8 days, InvivoGen's Anti-mCTLA4-mlgG2a or Competitor X's Anti-mCTLA4-mlgG2b mAb was administered intraperitoneally (IP) into the mice in sterile PBS (200 µg/mouse). Following this, IP injections were performed twice a week for 3 weeks (red arrows). Tumor growth was monitored for 60 days.

PRODUCT	DESCRIPTION	QTY	CAT. CODE
<i>In vivo</i> depletion mAbs			
Anti-mCD3-mlgG2a InvivoFit™ NEW	145-2C11-derived mouse mAb against murine CD3	1mg / 10 mg	mcd3-mab10-1 / -10
Anti-mCD4-mlgG2a InvivoFit™ NEW	GK1.5-derived mouse mAb against murine CD4	1mg / 10 mg	mcd4-mab10-1 / -10
Anti-mCD8-mlgG2a InvivoFit™ NEW	YTS169.4-derived mouse mAb against murine CD8	1mg / 10 mg	mcd8-mab10-1 / -10
<i>In vivo</i> immune checkpoint blocking mAbs			
Anti-mCTLA4-mlgG2a InvivoFit™	9D9-derived mouse mAb against murine CTLA-4	1mg / 10 mg	mctla4-mab10-1 / -10
Anti-mPD-1-mlgG1e3 InvivoFit™	RMP1-14-derived mouse mAb against murine PD-1	1mg / 10 mg	mpd1-mab15-1 / -10
Anti-PD-L1-mlgG1e3 InvivoFit™	Murinized atezolizumab mouse mAb against PD-L1	1mg / 10 mg	pdl1-mab15-1 / -10
<i>In vivo</i> tumor-associated antigen mAbs			
Anti-mgp75-mlgG2a InvivoFit™	TA99-derived mouse mAb against murine gp75	1mg / 10 mg	mgp75-mab10-1 / -10
<i>In vivo</i> isotype controls			
Anti-β-Gal-mlgG1e3 InvivoFit™	Mouse IgG1e3 isotype control	1mg / 10 mg	bgal-mab15-1 / -10
Anti-β-Gal-mlgG2a InvivoFit™	Mouse IgG2a isotype control	1mg / 10 mg	bgal-mab10-1 / -10

InvivoGen offers **bulk quantities** of mAbs in our mouse anti-mouse collection **upon request**. Please [contact us](#) for more information.

To learn more about our mouse anti-mouse mAb collection visit:

www.invivogen.com/mouse-anti-mouse-mabs

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
Fig 1: Infection by SARS-CoV-2 Spike pseudotyped lentiviral particles. HEK-Blue™ hACE2, A549-hACE2, and A549-hACE2-TMPRSS2 cells were cultured in the presence of SARS-CoV-2 Spike-D614 (HEK293) or Spike-G614 (A549) pseudotyped GFP lentiviral particles. After 72 hours, the transduction efficiency of the Spike-pseudotyped GFP particles was evaluated by fluorescence microscopy.



Fig 2: Neutralization ability of Anti-SARS-CoV-2 Spike-RBD mAbs. SARS-CoV-2 spike-RBD targeting mAbs (clones H4 & B38) were incubated with Spike (G614)-pseudotyped lentiviral particles for 1 hour at 37°C. Following this, HEK-Blue™ hACE2 cells were added and incubated for 72 hours. Infection (GFP fluorescence) was then measured using flow cytometry. Data are presented as % inhibition of infection compared to a no antibody control.

Application of cell-based assays for SARS-CoV-2 research

PRODUCT	QTY	CAT. CODE
HEK-Blue™ hACE2 Cells	3-7 x 10 ⁶ cells	hkb-hace2
A549-hACE2-TMPRSS2 Cells	3-7 x 10 ⁶ cells	a549-hace2tpsa
A549-hACE2 Cells (control)	3-7 x 10 ⁶ cells	a549-hace2

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For a complete list of SARS-CoV-2 related research tools visit:

www.invivogen.com/covid-19

InvivoGen

Europe Tel: +33 562 71 69 39
USA Tel: +1 888 457 5873
Asia Tel: +852 3622 3480

Fax: +33 562 71 69 30
Fax: +1 858 457 5843
Fax: +852 3622 3483

info.eu@invivogen.com
info@invivogen.com
info.hk@invivogen.com