
# INFLAMMASOMES

InvivoGen *infocus*

PRACTICAL GUIDE

## UNDERSTANDING INFLAMMASOMES

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# InvivoGen *infocus*

## A PRACTICAL GUIDE TO UNDERSTANDING INFLAMMASOMES

4

### **INFLAMMASOME ACTORS**

Inflammasome sensors  
Adaptors  
Effectors

6

### **INFLAMMASOME CELLULAR ASSAYS**

Monitoring NLRP3-dependent responses  
Monitoring NLRC4-dependent responses  
Monitoring ASC-dependent responses  
Monitoring CASP-1-dependent responses  
Monitoring CASP-4-dependent responses

8

### **INFLAMMASOME ACTIVATION AND ASSEMBLY**

A "two-step" activation model  
NLRP3 inflammasome  
NLRC4/NAIP inflammasome  
CASP-11 & CASP-4/5 non-canonical inflammasomes

10

### **INFLAMMASOME ACTIVATION & RESPONSES AT A GLANCE**

Poster

12

### **INFLAMMASOME INDUCERS**

Canonical inflammasome inducers  
Non-canonical inflammasome inducers

14

### **PYROPTOSIS AND GASDERMIN D**

#### **GASDERMIN D CELLULAR ASSAY**

Monitoring GSDMD-dependent responses

15

### **INFLAMMASOME REGULATION**

#### **PRIMING AGENTS**

LPS  
Poly(I:C)  
Pam3CSK4

16

### **ASSESSING INFLAMMASOME ACTIVATION**

Monitoring ASC-specks in real-time  
Monitoring pyroptotic cell death  
Monitoring IL-1 $\beta$  and IL-18 secretion

18

### **INFLAMMASOME INHIBITORS**

Caspase inhibitors  
Multi-target inhibitors  
NLRP3 inhibitors

19

### **INFLAMMASOMES IN DISEASE REFERENCES**

## ABBREVIATIONS

**AIM2**: Absent-in-melanoma-2

**ATP**: Adenosine tri-phosphate

**ASC**: Apoptosis-associated speck-like protein

**CARD**: Caspase recruitment domain

**CASP**: Caspase

**DAMPs**: Danger-associated molecular patterns

**GBPs**: Guanylate binding proteins

**GSDM**: Gasdermin

**HMGB1**: High mobility group B1 protein

**IFN**: Interferon

**IL-1**: Interleukin 1

**IRF**: Interferon regulatory factor

**IRGB10**: Interferon response gene B10 protein

**LDH**: Lactate dehydrogenase

**LRR**: Leucine-rich repeat

**LPS**: Lipopolysaccharide

**NAIP**: NLR family apoptosis inhibitory protein

**NEK7**: NimA-related protein kinase 7

**NF- $\kappa$ B**: Nuclear factor kappa-light-chain enhancer of activated B cells

**NBD**: Nucleotide-binding domain

**NLR**: Nucleotide-binding domain and leucine-rich repeat

**NLRP3**: NLR family PYD domain containing protein 3

**NLRC4**: NLR family CARD domain containing protein 4

**NOD**: Nucleotide-binding oligomerization domain

**OMVs**: Outer membrane vesicles

**PAMPs**: Pathogen-associated molecular patterns

**PRRs**: Pattern recognition receptors



**PTM**: Post-translational modification

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## A BRIEF HISTORY OF INFLAMMASOME BREAKTHROUGHS

**6 (AT LEAST) INFLAMMASOME SENSORS**


- NLRP1
- NLRP3
- NLRP4
- PYRIN
- CASPASE II-4/5
- AIM2

**3 PROTEINS IN CANONICAL ASSEMBLY**

SENSOR + ADAPTOR + CASPASE-1

**2 SIGNALS FOR INFLAMMASOME FORMATION**

PRIMING + ACTIVATION

**70**

**STUDIES IN** [clinicaltrials.gov](https://clinicaltrials.gov)

**InvivoGen**

**OFFERS >60**

CELLULAR ASSAYS, REAGENTS, AND GENE TOOLS TO ACCOMPANY YOUR SUCCESSFUL RESEARCH

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
## Monitoring NLRP3-dependent responses

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#### Functional characterization of THP1-KO-NLRP3 Cells.

THP1-KO-NLRP3 cells and their parental cell line, THP1-Null2, were primed with LPS-EK (1 μg/ml) then stimulated with Nigericin (5 μM), Alum Hydroxide (150 μg/ml), transfected Poly(dA:dT) (1 μg/ml), or E. coli OMVs (100 μg/ml). After overnight incubation, IL-1β secretion was assessed using HEK-Blue™ IL-1β sensor cells and the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 630 nm.

PRODUCT	QTY	CAT. CODE
THP1-KO-NLRP3 Cells	3-7 x 10 <sup>6</sup> cells	thp-konlrp3z
THP1-Null2 Cells	3-7 x 10 <sup>6</sup> cells	thp-nullz

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
## Monitoring NLRC4-dependent responses

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#### Functional characterization of THP1-KO-NLRC4 and RAW-ASC KO-NLRC4 Cells.

(A) THP1-KO-NLRC4 cells and their parental cell line, THP1-Null2, were primed with LPS-EK (1 μg/ml). (B) RAW-ASC KO-NLRC4 cells and their parental cell line RAW-ASC, were primed with Pam3CSK4 (100 ng/ml). (A and B) After priming, the cells were stimulated with Poly(dA:dT) (1 μg/ml), Needle-Tox (4 ng/ml) or Rod-Tox (2 μg/ml). After overnight incubation, the secretion of human or mouse IL-1β was assessed using (A) the HEK-Blue™ IL-1β cell based assay or (B) an ELISA assay, respectively.

PRODUCT	QTY	CAT. CODE
THP1-KO-NLRC4 Cells	3-7 x 10 <sup>6</sup> cells	thp-konlrc4z
RAW-ASC KO-NLRC4 Cells	3-7 x 10 <sup>6</sup> cells	raw-konlrc4
RAW-ASC Cells	3-7 x 10 <sup>6</sup> cells	raw-asc


Sharma B. et al. 2019. Methods Mol. Biol. 1682:135-47. Assessing NLRP3 inflammasome activation by nanoparticles.

Jäger E. et al. 2020. Nat. Commun. 11(1):4243. Calcium-sensing receptor-mediated NLRP3 inflammasome response to calciprotein particles drives inflammation in rheumatoid arthritis.

## Monitoring ASC-dependent responses

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#### Functional characterization of THP1-KO-ASC Cells.


THP1-KO-ASC cells and their parental cell line, THP1-Null2, were primed with LPS-EK (1 µg/ml) and then stimulated with Nigericin (5 µM), Alum Hydroxide (150 µg/ml), transfected Poly(dA:dT) (1 µg/ml), or E. coli OMVs (100 µg/ml). After 24h activation, IL-1β secretion was assessed in the culture supernatant using HEK-Blue™ IL-1β sensor cells and the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 630 nm.

PRODUCT	QTY	CAT. CODE
THP1-KO-ASC Cells	3-7 x 10 <sup>6</sup> cells	thp-koascz

## Monitoring CASP-1-dependent responses

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#### Functional characterization of THP1-defCASP1 Cells.

THP1-defCASP1 cells and their parental cell line, THP1-Null, were primed with LPS-EK (1 µg/ml) prior to stimulation with Nigericin (5 µM), MSU crystals (250 µg/ml), or E. coli OMVs (100 µg/ml). After 24h activation, IL-1β secretion was assessed in the culture supernatant using HEK-Blue™ IL-1β sensor cells and the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 630 nm.

PRODUCT	QTY	CAT. CODE
THP1-defCASP1 Cells	3-7 x 10 <sup>6</sup> cells	thp-dcasp1
THP1-Null Cells	3-7 x 10 <sup>6</sup> cells	thp-null


## STUDYING SPECIFIC INFLAMMASOME ACTORS

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## Monitoring CASP-4-dependent responses

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#### Functional characterization of THP1-KO-CASP4 Cells.

THP1-KO-CASP4 cells and their parental cell line, THP1-Null2, were primed with LPS-EK (1 µg/ml) then stimulated with Nigericin (5 µM; canonical inducer) or transfected LPS-EK (5 µg/ml, non-canonical inducer). After 6h incubation, IL-1β secretion was assessed using HEK-Blue™ IL-1β sensor cells and the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 630 nm.

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# INFLAMMASOME ACTIVATION AND ASSEMBLY

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
# INFLAMMASOME ACTIVATION & RESPONSES AT A GLANCE


**PRE-PRIMING**  
(for non-canonical inflammasomes)

**PRIMING (signal 1)**









## ACTIVATION (signal 2)

## OUTCOMES

## OVERVIEW OF INFLAMMASOME SENSORS AND ACTIVATORS

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

#### Induction of IL-1 $\beta$ secretion by monocytes upon treatment with Nigericin.

Human THP-1 monocytes were primed with LPS-EK (1  $\mu$ g/ml) prior to incubation with increasing concentrations of Nigericin. After 6h and 24h activation, IL-1 $\beta$  secretion was assessed in the culture supernatant using HEK-Blue™ IL-1 $\beta$  sensor cells and the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 630nm.

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#### Induction of IL-1 $\beta$ secretion by monocytes upon treatment with MSU crystals.


Human THP-1 monocytes were primed with LPS-EK (1  $\mu$ g/ml) prior to incubation with increasing concentrations of MSU crystals. After 6h and 24h activation, IL-1 $\beta$  secretion was assessed in the culture supernatant using HEK-Blue™ IL-1 $\beta$  sensor cells and the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 630nm.

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#### Induction of IL-1 $\beta$ secretion by monocytes upon transfection with Poly(dA:dT).

Human THP-1 monocytes were primed with LPS-EK (1  $\mu$ g/ml) prior to transfection with increasing concentrations of Poly(dA:dT). After 6h and 24h activation, IL-1 $\beta$  secretion was assessed in the culture supernatant using HEK-Blue™ IL-1 $\beta$  sensor cells and the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 630nm.

## TIPS

### INDUCING INFLAMMASOME ACTIVATION

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
## Non-canonical inflammasome inducers

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#### Induction of IL-1β secretion by monocytes upon transfection with LPS-EB.

Human THP-1 monocytes were primed with LPS-EK (1 μg/ml) prior to transfection with increasing concentrations of LPS-EB. After 6h and 24h activation, IL-1β secretion was assessed in the culture supernatant using HEK-Blue™ IL-1β sensor cells and the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 630 nm.

Tweedell R.E. *et al.*, 2020. Nat. Protocols. 15(10):3284-3333. A comprehensive guide to studying inflammasome activation and cell death.



Jain A. *et al.*, 2020. Nat. Immunol. 21:65-74. T cells instruct myeloid cells to produce inflammasome-independent IL-1β and cause autoimmunity.

Zheng M. *et al.*, 2020. Cell. 181:674-687.e13. Caspase-6 Is a Key Regulator of Innate Immunity, Inflammasome Activation, and Host Defense.

Samir P. *et al.*, 2019. Nature. 573:590-594. DDX3X acts as a live-or-die checkpoint in stressed cells by regulating NLRP3 inflammasome.

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#### (A) Functional characterization of E. coli OMVs.

THP1-HMGB1-Lucia™ cells (see page 16) were incubated with increasing concentrations of *E. coli* OMVs. After 24h activation, secretion of HMGB1 was assessed by measuring the Lucia luciferase activity in the culture supernatant using QUANTI-Luc™ detection reagent.

(B) Image of InvivoGen's *E. coli* OMVs by transmission electron microscopy (80 kV).

## KEY POINT

### DIFFERENTIAL CASPASE 11/4/5 ACTIVATION BY VARIOUS LPS STRUCTURES

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
## GASDERMIN D CELLULAR ASSAY

### Monitoring GSDMD-dependent responses

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**Functional characterization of THP1-KO-GSDMD and RAW-ASC KO-GSDMD Cells.** THP1-Null2 parental and THP1-KO-GSDMD cells (A) or RAW-ASC parental and RAW-ASC KO-GSDMD cells (B) were primed with LPS-EK or Pam3CSK4, prior to stimulation with Nigericin (5 μM) or E. coli outer membrane vesicles (OMVs) (100 μg/ml). After 6h activation, human IL-1β secretion was assessed in the culture supernatant using HEK-Blue™ IL-1β sensor cells (see page 7) (A). Murine IL-1β secretion was assessed in the culture supernatant using an ELISA assay (B).

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


# ASSESSING INFLAMMASOME ACTIVATION

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Visualization of ASC speck formation by fluorescence microscopy.

A. Procedure for inducing ASC::GFP expression in THP1-ASC-GFP cells. B. ASC speck visualization in THP1-ASC-GFP cells after priming with 1  $\mu$ g/ml LPS-EK for 3 h, and activation with 250 ng/ml Poly(dA:dT) for 1 to 3 h. In most cells, only one speck forms upon inflammasome activation (arrows). Scale bar: 50  $\mu$ m.

PRODUCT	QTY	CAT. CODE
THP1-ASC-GFP Cells	3-7 x 10 <sup>6</sup> cells	thp-ascgfp

## TIPS

### EVALUATION OF INFLAMMASOME ACTIVATION

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


Luminescent quantification of inflammasome-induced pyroptosis.

A. Assay principle for THP1-HMGB1-Lucia™ cells. B. THP1-HMGB1-Lucia™ cells were primed with LPS-EK (1  $\mu$ g/ml) for 3 h and then incubated with the inflammasome inducer Poly(dA:dT) (0.5  $\mu$ g/ml). Lucia luciferase activity (pink curve, left axis) and LDH release (blue curve, right axis) in the supernatant were quantified at 2, 3, 4, 5 and 6 h post-induction.

PRODUCT	QTY	CAT. CODE
THP1-HMGB1-Lucia™ Cells	3-7 x 10 <sup>6</sup> cells	thp-gb1lc

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
**Dose response of HEK-Blue™ IL-1 $\beta$  and HEK-Blue™ IL-1R cells to human and murine IL-1.** (A) HEK-Blue™ IL-1 $\beta$  and (B) HEK-Blue™ IL-1R cells were incubated with increasing concentrations of recombinant hIL-1 $\alpha$ , hIL-1 $\beta$ , mL-1 $\alpha$ , or mL-1 $\beta$ . After overnight incubation, SEAP activity was measured in the supernatant using QUANTI-Blue™ Solution.


## ADVANTAGES

### USING INVIVOGEN'S CYTOKINE SENSOR CELLS

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**Dose response of HEK-Blue™ IL-18 cells to human and murine IL-18.** HEK-Blue™ IL-18 cells were incubated with increasing concentrations of recombinant human or murine IL-18. After overnight incubation, levels of NF- $\kappa$ B/AP1-induced SEAP in the supernatant were determined using the SEAP detection reagent QUANTI-Blue™ Solution. The optical density (OD) was read at 655 nm.

PRODUCT	QTY	CAT. CODE
HEK-Blue™ IL-1 $\beta$ Cells	3-7 x 10 <sup>6</sup> cells	hkb-il1bv2
HEK-Blue™ IL-1R Cells	3-7 x 10 <sup>6</sup> cells	hkb-il1r
HEK-Blue™ IL-18 Cells	3-7 x 10 <sup>6</sup> cells	hkb-hmil18

Martine P. et al. 2019. Cell Death Dis. doi: 10.1038/s41419-019-1491-7. HSP70 is a negative regulator of NLRP3 inflammasome activation.

Metho S. et al. 2019. Mol. Cell. doi: 10.1016/j.molcel.2018.11.08. The Crohn's disease risk factor IRGM limits NLRP3 inflammasome activation by impeding its assembly and by mediating its selective autophagy.

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They trust **InvivoGen**

Hafner-Bratkovič I. et al., 2018 Nat. Commun. 59:5182. NLRP3 lacking the leucine-rich repeat domain can be fully activated via the canonical inflammasome pathway.  
Irmischer S. et al., 2019 Nat. Commun. 10:2961. Serum FHR1 binding to necrotic-type cells activates monocyctic inflammasome and marks necrotic sites in vasculopathies.  
da Costa L.S. et al., 2019 Cell Death Dis. 10:346. RNA viruses promote activation of the NLRP3 inflammasome through cytopathogenic effect-induced potassium efflux.

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