



ISA Plate Seal / FSA Plate Seal



Manufactured by LINKSu

96-well plate seal; ISA plate seal

Unique self-closing idea prevents solvent evaporation



Data: Comparison of water and acetonitrile



150 µl of water or acetonitrile in 96-well plate

Compatible with below auto-sampler systems

- Nano space by Shiseido
- Prominence by Shimadzu
- ACQUITY UPLC by Waters
- PAL-HTC by CTC
- 1100, 1200 and 1290 series by Agilent

Product specification

Material: Silicon, PET
 Operating temperature: -80 to 100 °C
 Size: 80 mm x 122 mm


Ordering information

Product name	Storage	Product number	PKG Size
ISA Plate Seal	RT	LKS-BCS96-80122	100 sheets


96-well plate seal; FSA plate seal

High resistance to organic solvents eliminates elution of adhesive.

Unique embossed surface prevents needle clogging.



Data: Highly resistant to organic solvents such as DMSO, acetonitrile or methanol.



Compatible with below auto-sampler systems
 ■ 6400 Series Triple Quadrupole LC/MS System by Agilent
 ■ Prominence SIL-20A by Shimadzu
 ■ Nexera X2 SIL-30AC by Shimadzu

Product specification
 Material: Polyolefin, PET
 Operating temperature: -80 to 80 °C
 Size: 80 mm x 122 mm

Ordering information

Product name	Storage	Product number	PKG Size
FSA plate seal	RT	LKS-BC001	100 sheets

For research use only, not intended for diagnostic or drug use.

NACALAI TESQUE, INC.

Nijo Karasuma, Nakagyo-ku, Kyoto 604-0855, JAPAN
 TEL : +81(0)75-251-1730
 FAX : +81(0)75-251-1763
 Website : <http://www.nacalai.com>
 E-mail : info.intl@nacalai.com